

# centrix<sup>®</sup>


MAKING DENTISTRY EASIER<sup>®</sup>

PERFECT  
PREVENTION

EFFICIENT  
DELIVERY

ACCURATE  
IMPRESSIONS

## BEAUTIFUL RESULTS!



Celebrating

50  
Years

“Like most practicing dentists, I experienced problems with products that made the delivery of quality dentistry more difficult than it needed to be. My first products, AccuDose® tubes and the C-R® syringes, changed the use of dental materials and made dentistry easier. 50 years later, I am proud that Centrix continues to introduce products and ideas that fulfill our mission.”

**Making Dentistry Easier.**

~ Dr. William Dragan



**MAKING DENTISTRY EASIER!**

Today, Centrix manufactures and engineers hundreds of dental products focused primarily on general restorative and preventive dentistry, the heart and soul of most dental practices. As a leader in materials and delivery systems, Centrix continues to introduce inventive dental systems certain to become the new standard in dentistry.

**OUR FIRST INNOVATIONS!**



**Mark IIIp™ Gun & Snap-Fit™ Gun**

Easy, direct-to-prep delivery of today's dental materials. Made in the USA and Used in every dental office!



**Benda® Brush Applicators**

One-piece, bendable head, precision brushes. Made in the USA. Can you imagine dentistry without them?



**AccuDose® & NeedleTube®**

Direct-to-prep done better! Dispense exactly the material needed with no waste or cross-contamination. Centrix Tubes & Plugs are made in the USA and truly changed the industry!

**OUR NEWEST INNOVATIONS!**



**NoCord™ VPS**

Industry-First: Hemostasis and retraction built-in! Award-winning, one-step, self-retracting, impression material. Better marginal detail and better fluids control.



**Encore® D/C + Connexio®**

Making post-end treatment easier and more assured. For restorations that last! Increased bond strengths reduce polymerization shrinkage, virtually eliminating micro-leakage.



## PREVENTION FOR LIFE®

Caries rates in kids are decreasing! However, caries rates in adults are increasing, rapidly, as patients decline to pay for treatment not covered by insurance – or are not offered preventive care, as a matter of course.

The Centrix Prevention for Life program provides tools and information to help offices update their preventive treatment protocols and how they discuss them with patients.

**Webinars and in-office Lunch & Learn** programs that provide CE credits.

**Risk Assessment Tools** that help you quickly, easily and visually share with your patients your diagnosis. **Prescription and Treatment Program** that increases patient health (as well as practice health). **Products** that are easy to place and your patients will love!

**At Centrix, we believe that structured Risk Assessment quickly results in diagnosis of risk and a prescription for that patient (rinse, varnish, etc.), treatment plan and frequency for each patient, and easy patient acceptance!**

Start a Prevention For Life® program today! Visit [cariesprevention.com](http://cariesprevention.com) or call 800.235.5862



# ORAL CARE AND PREVENTION

Centrix makes prevention easier for you and your patients



## FluoroDose®

Award-Winning Fluoride Varnish with Xylitol



- Tastes Better! Feels Better! Better Priced!
- Application takes less than a minute.
- Unit-dose delivery means easier use, easier cleanup.

### TOP FLUORIDE VARNISH, 6 YEARS IN A ROW!



### DENTAL ADVISOR TOP AWARD: FLUORIDE VARNISH



## LoliCaine®

Topical Anesthetic | 20% Benzocaine Gel

- 4 Great flavors: Bubble Gum, Cherry, Mint, Piña Colada.
- Rapid onset with no alcohol aftertaste!
- Unit-dose LolliPack® saves time and reduces cross-contamination.
- The LolliPack is the perfect “take-home pack” for post-op discomfort.



## E-Lieve™

Vitamin E Treatment | Made in the USA

E-Lieve helps to relieve:

- Chapped lips before, during and after dental procedures.
- Dryness and irritation when using rubber dams.
- Sore tissue from deep scaling procedures.
- Discomfort from wearing new dentures.
- Irritation caused by orthodontic brackets.

**A Faster, Easier Solution for Moisturizing and Soothing Minor Soft Tissue Irritation.**



## Champ®

Pit & Fissure Sealant | Glass-Filled Acrylate Resin

- **Bioactive:** hydrophilic resin forms a strong resin-hydroxyapatite complex that seals out micro-leakage.
- **Moisture-friendly formula:** tolerates moisture making it easier to bond to the tooth.
- **Margin Integrity:** lightly filled (37%) resin is durable and long-lasting; adapts to and seals the margins, helping to eliminate micro-leakage.



## Exposé®

Caries Indicator on a Stick

- Easily identify and remove all infected dentin.
- Single-patient applicators make tray set-up easy.
- No syringe, no sheath, no tips to replace.
- No cross-contamination, no messy spills, no stains!
- Made in the USA.

# IMPRESSIONING

Versatile choices for easier, safer impressions

UNIQUE AND COST-EFFECTIVE SOLUTIONS FOR ALL IMPRESSION TECHNIQUES

DENTISTRY TODAY  
**TOP Innovative Products 2019**

Dental PRODUCT SHOPPER  
**BEST PRODUCT**



## NoCord™ VPS

One-Step, Self-Retracting Impressioning System

- The safest, most unique VPS available today.
- Built-in hemostasis – wash stops bleeding and fluid flow for clear margins, and provides exceptional detail.
- MegaBody™ Tray Material creates sulcular expansion. Enables complete margin capture, usually without retraction cord. Adds extra hemostasis when you do use cord.



Made in the USA



## Access GOLD®

- Low-viscosity wash material, syringes easily into sulcus and around prep.
- High-contrast color makes it easy to verify all margins have been captured.

## Access HEAVY®

- High-viscosity tray material, and the perfect complement to Access GOLD.
- Easy to place, easy to remove.

## Access MONOPHASE®

- Non-slumping, medium-body material.
- Perfect for edentulous impressions or as a tray material.

## Access PUTTY®

- Non-oily and easy to mix.
- For use as a tray material or for easy fabrication of stents, matrices and more.

EASIER, FASTER BITE REGISTRATION

Dental PRODUCT SHOPPER  
**BEST PRODUCT**



## Access BLUE®

VPS Quick Bite Registration Material

- Mousse-like viscosity for easy, fast application.
- Sets hard, like laboratory plaster in just 30 seconds.
- Most accurate bite registration for intercuspal and centric relationships.
- Cut and trim the impression without it crumbling.
- Made in the USA.



## 3.n.One™

Impression Trays For Dual-Arch Technique

- Unique undercut ridges lock impression materials in place.
- High rigidity with non-tear, flexible webbing ensures maximum accuracy.
- Three different shapes to meet all applications.

BETTER RETRACTION WHEN YOU NEED IT

REALITY Five Star Award  
★★★★★  
REALITY'S CHOICES



## Access® Edge

High-Viscosity Kaolin Clay with 15% Aluminum Chloride

- A thicker, more viscous paste that provides more 'push' of the soft tissue.
- Uses standard composite syringe; no glass vials.
- Dual action – stops bleeding, and absorbs moisture.
- Made in the USA.
- GingiCap™ Compression Caps included.

Two Kaolin Clay formulations to meet every need!



## AccessFLO®

Low Viscosity Kaolin Clay with 15% Aluminum Chloride

- Flowable paste that more easily penetrates the sulcus.
- Dual action – 15% aluminum chloride plus kaolin clay stops bleeding and absorbs moisture.
- Made in the USA.
- GingiCap™ Compression Caps included.

THE DENTAL ADVISOR  
+++++½



## GingiTrac™

VPS Gingival Retraction System

- Makes digital dentistry easier and better.
- Unique VPS material gently retracts tissue yet removes cleanly with no debris or particles that can distort digital images.
- Ammonium aluminum sulfate to control bleeding and oozing.
- Made in the USA.
- GingiCap™ Compression Caps included.

# POST-ENDO TREATMENT

## Centrix products provide a superior approach

### ROOT CANAL TREATMENT



#### Calciject® Calcium Hydroxide Paste for Root Canal Therapy

- No messy mixing of powder and liquid.
- Direct injection from package to mouth.
- Single-dose tips for maximum asepsis.
- Bendable 20 ga NeedleTube® tips for easy canal access.
- Superior results compared to spiral paste fillers.
- Bacteriostatic and radiopaque.
- Made in the USA.



#### Connexio® Dual-Cure, Self-Etch Bond for Core Buildups

- Designed for the demands of deep core build-ups where light cannot reach.
- Dual-cure adhesive system provides the highest possible bond strength.\*
- Nanoparticle reinforced.
- Water solvent (not acetone or ethyl alcohol); makes Connexio more tolerant of humid canal conditions than other adhesives.

\*Migliau, G; Piccoli, L; Di Carto, et al, Comparison between three glass fiber post cementation techniques. (Ann. Stomatol (Roma)); 2017 Jan-Mar; 8(1): 29-33. Published online 2017 Jul 3 doi [10.11138/ads/2017.8.1.029].

More dentists choose **Encore D/C** and **Connexio** for a bond that can't be broken!



#### Encore® D/C MiniMix® Dual-Cure, Core Buildup Composite

- Non-slumping formula stays where it's placed, will not run and does not require a matrix or core form.
- Exclusive MiniMix® cartridge (0.5 mL) controls costs while eliminating waste and cross-contamination.
- Cuts like dentin to reduce ditching.
- NeedleTube® mixing tip delivers material to bottom of the post hole.
- Made in the USA.
- **Also available in traditional 25 mL automix cartridges.**



### ESSENTIAL TOOLS FOR STRUCTURAL RESTORATIVE INTEGRITY




#### Absolute™ Universal Resin Cement

- Self-etch, self-adhesive, dual-cure for faster, easier placement. No primer or conditioner needed.
- Definitive cementation that lasts.
- Designed especially for feldspathic, high leucite, lithium disilicate, alumina and zirconia materials.
- Made in the USA.

#### 3 Shades = A natural-looking restoration

- **Translucent:** Provides a natural esthetic
- **Opaque (White):** Completely blocks underlying colors/show-through
- **Universal (Opaque Tooth Shade):** Completely blocks underlying colors/show-through



#### Onyx™ 40% Phosphoric Acid Etching Gel

- Highest strength, most effective etchant.
- Black silicon carbide particles provide superior contrast.
- Better control and easier placement.
- Stays where placed, will not run.
- Particles can be agitated on tooth to mechanically scrub and clean the surface.
- No dye, does not stain, rinses completely.
- Made in the USA.

# TEMPORIZATION

Complete solutions designed to work better together

THE RIGHT FILLING MATERIAL FOR EVERY CLINICAL SITUATION



Short-Term Temps

## Tempit® Moisture-Activated Temporary Filling & Sealing Material

- Ideal for endo access openings and short-term temporaries.
- Moisture-activated, sets in 5 minutes or less.
- Expands slightly when set, seals out microleakage and bacterial ingress.
- Contains calcium sulfate and zinc oxide.
- Significant bacteriocidal activity up to 14 days.<sup>1</sup>

<sup>1</sup>Antibacterial Properties of Temporary Filling Materials, JOE, Volume 32, Number 3, March 2006



Short-Term Temps with Eugenol

## Tempit®-E Moisture-Activated Temporary Filling & Sealing Material with Eugenol

- Palliative eugenol formulation for temporary filling use, especially in cases of reversible pulpal inflammation.
- Anti-inflammatory and local anesthetic effects helps facilitate pulpal healing.
- Moisture activated, sets in five minutes.
- Expands slightly when set to seal out microleakage and bacterial ingress.
- Can be used as a liner in deep preparations.



Indirect Temp Restorations

## Tempit® L/C Light-Activated Temporary Material for Inlays & Onlays

- Perfect for indirect restorations, inlays/onlays.
- Seals implant screw access openings.
- Sets firm yet flexible.
- Easily removed with an explorer without damaging the margins.
- Microfilled polymer resin.
- Light-activated, sets in 40 seconds or less.
- No etching or bonding agent required.



Longer-Term Temps

## Tempit® Ultra-F Light-Activated, Longer-Term Temporary Material

- Perfect for longer term temporaries.
- Rigid, harder wearing material, with fluoride release.
- High compressive strength for bruxism patients.
- Ideal for tooth-colored restorations in deciduous teeth.
- Less soluble than reinforced ZOE intermediate restorative material.
- Light-activated, sets in 40 seconds or less.
- No etching or bonding agent required.
- **Tempit Ultra: without fluoride**



## TWO BETTER TEMPORARY CEMENTS



## Access® Automix Temporary Cement

- Perfect for longer-term temporary use.
- Low film thickness.
- Flows and places easily.
- Semi-elastic set.
- Safe, easy adhesive.
- Calcium hydroxide provides bacteriostatic effect on mold.
- Eugenol-free, won't interfere with bonding agents or other resin materials.



Actual size.

## NoMIX® Short-term, Temporary Cement

- Easy to use, dispense, clean and remove from prep.
- Low film thickness, harder set.
- Tube applicators are ideal for both chair-side application and for patient take-home.
- Expands slightly upon setting for a tight margin seal.
- Moisture activated, no need to dry.
- **Also available in syringe delivery for chair-side application.**



## SUPERIOR BUILDUP MATERIAL





## Access Crown® Temporary Crown & Bridge Material

- **More:** shades (7), speed, strength and savings.
- **Ideal handling characteristics:** ready to remove from mouth in 60 seconds.
- Low shrinkage, low heat generation.
- High compressive, flexural strength.
- Highly aesthetic – polishes easily, beautifully.



Natural Esthetics, High Fluorescence, Polishes Easily.



# DELIVERY SYSTEMS

Control costs and eliminate cross-contamination

APPLY MATERIALS EASILY AND MORE PRECISELY

## Benda® Brush

The Original Precision Disposable Brush Applicator

- When you need to apply more than a drop.
- Bend head to any angle with no bounce-back!
- Unique non absorbing fibers have a greater carrying capacity, ensuring better spreading.
- Regular point for large applications.
- Mini point for smaller spaces.
- Made in the USA.



Benda Twin with snap-apart handle. The lowest cost alternative!

More sizes, more styles and more colors means the right applicator for every material and every need!

## Benda® Micro

Precision Disposable Fiber Applicator

- Perfect for when you need only a drop.
- Regular point (2.0 mm).
- Fine point (1.5 mm).
- Extra-Fine point ( 1.0 mm).
- Non-linting design will not leave fibers behind in restorations or sulcuses.
- Made in the USA.

## Benda® Twin

Double-ended, Twin-style Applicator

- When you are using more than one material.
- Ideal for 2-step procedures.
- Extended reach.
- Made in the USA.

Benda Brush

Benda Micro



DELIVER MATERIALS MORE EFFECTIVELY



More Usable Material Per Cartridge!

■ Mixer Waste (g)  
■ Usable Material (g)



## SuperMixer™ AutoMix Mixing Tips

- High-efficiency mixer, up to 38% reduced material waste.
- Up to 50% more units from every cartridge.
- More consistent mixing of any material.
- Color-coded hub for easy replacement of old-style tips.
- Replaces any 1:1 mixing tip – same color for same color.

## AccuDose® & NeedleTube® Delivery Systems

- Dispense exactly the material needed.
- Prevent cross-contamination.
- Ready to use, dispense directly into the prep.
- No waste, no contamination.
- No more 'composite keeper'
- Compatible with all C-R Syringes and guns.
- Made in the USA.

Less waste in every tip!

PREVENTION

**FLUORODOSE® UNIT DOSE, 120-CT.**

- 360086 Mint, 0.3 mL LolliTray™
- 360087 Bubble Gum, 0.3 mL LolliTray™
- 360107 Melon, 0.3 mL LolliTray™
- 360110 Cherry, 0.3 mL LolliTray™
- 360137 Caramel, 0.3 mL LolliTray™

**CHAMP®**

**Unit Dose**

- 360116 16-ct., 0.25 mL LolliPack
- 360120 64-ct., 0.25 mL LolliPack
- 310033 Champ Sealant Syringe Kit 4x1.2 mL, 20 tips
- 310034 Champ Sealant Syringe Single 1x1.2 mL, 5 tips
- 310032 22 ga Luer Lock Needle Tips, 20-ct.
- 640001 22 ga Luer Lock Needle Tips, 100-ct.

**LOLLICAINE® UNIT DOSE, 120-CT.**

- 360088 Bubble Gum, 0.3 mL LolliPack
- 360090 Cherry, 0.3 mL LolliPack
- 360092 Mint, 0.3 mL LolliPack
- 360094 Piña Colada, 0.3 mL LolliPack

**E-LIEVE™**

- 310171 50 x 0.5 mL Applicators

**EXPOSÉ®**

- 310334 Fine, 100-ct., Blue
- 310330 Fine, 576-ct., Blue
- 310335 Extra Fine, 100-ct., Blue
- 310331 Extra Fine, 576-ct., Blue

IMPRESSIONING

**NOCORD™ VPS**

- 310280 VPS Intro Trial Kit
- 310285 VPS System Intro Kit
- 310282 Wash Material, 2-ct.
- 310284 Wash Material, 10-ct.
- 310281 MegaBody™ Tray Material, 2-ct.
- 310283 MegaBody Tray Material, 10-ct.
- 310286 MegaBody, 380 mL, 1-ct.
- 310287 MegaBody, 380 mL, 2-ct.

**ACCESS®**

- 360002 GOLD®, 2 x 50 mL
- 360016 HEAVY®, 2 x 50 mL
- 360017 MONOPHASE®, 2 x 50 mL
- 360015 PUTTY®, 2 x 280 mL Jars

**ACCESS RETRACTION PASTE**

- 360113 Access® Edge, 15-ct.
- 360114 Access® Edge, 60-ct.
- 640002 16 ga Luer Lock Needle Tip, 30-ct.
- 360122 AccessFLO®, 15-ct.
- 360123 AccessFLO®, 60-ct.

**GINGITRAC™**

- 360100 Refill Kit, 2 x 50 mL
- 360101 Standard Kit, 2 x 50 mL, 1:1/2:1 Gun
- 360102 Plus Kit, 4 x 50 mL
- 360068 MiniMix Starter Kit, 24-ct.
- 360069 MiniMix Refill Kit, 48-ct.

**ACCESS® BLUE**

- 360006 Bite Registration Kit, 2 x 50 mL

**3.N.ONE™ IMPRESSION TRAYS**

- 360011 Anterior, 35-ct.
- 360012 Posterior, 50-ct.
- 360013 Quadrant, 35-ct.
- 360014 Assorted, 40-ct.

**POST-ENDO**

**CALCIJECT®**

- 320003 Prefilled Tip, 0.20 g, 24-ct.

**CONNEXIO®**

- 310270 Self-Etch Bond, 3 mL Part A & 3 mL Part B

**ENCORE® D/C**

**MiniMix® Cartridge 0.5 mL**

- 310129 Natural, 12-ct.
- 310254 Natural, 36-ct.
- 310130 Contrast Blue, 12-ct.
- 310255 Contrast Blue, 36-ct.
- 700000 MiniMix Kit, 12-ct., Natural, with Snub-Nose Snap-Fit Syringe

**Automix Cartridge 25 mL, 10 mix tips**

- 310150 Natural, 1 x 25 mL
- 310152 Contrast Blue, 1 x 25 mL

**MiniMix® Nozzles**

- 910630 AutoMix 1:1 Mixing Nozzles, 10-ct.

**ABSOLUTE® UNIVERSAL RESIN CEMENT**

- 310113 Opaque, 1 x 5 g Syringe
- 310114 Translucent, 1 x 5 g Syringe
- 310015 Universal, 1 x 5 g Syringe
- 310016 Assorted, 3 x 5 g Syringe
- 210008 Blunt Mixing Tips, 50-ct.
- 210007 Tapered Mixing Tips, 50-ct.
- 210009 IntraOral Tips, 50-ct.

**ONYX™ ETCH**

- 310053 6 x 1.2 mL Syringes, 30 Needle Tips
- 310054 15 x 1.2 mL Syringe, 50 Needle Tips
- 310055 15 x 1.2 mL Syringe, 100 Needle Tips
- 310052 25 x 1.2 mL Syringe Only
- 310032 22 ga Luer Lock Needle Tips, 20-ct.
- 640001 22 ga Luer Lock Needle Tip, 100-ct.

TEMPORIZATION

**TEMPIT® LINE, PREFILLED TIPS**

- 310060 Tempit, 0.35 g, 30-ct.
- 310031 Tempit-E, 0.35 g, 30-ct.
- 310063 Tempit L/C, 0.25 g, 30-ct.
- 310064 Tempit Ultra-F, with fluoride, 0.20 g, 30-ct.
- 310065 Tempit Ultra, no fluoride, 0.20 g, 30-ct.
- 310067 Tempit Essentials Kit (5 ea. of Tempit, Tempit-E, Tempit L/C, Tempit Ultra-F)

**ACCESS® CROWN®, 76 g CARTRIDGE**

- 360029 Shade A3.5
- 360030 Shade A3
- 360031 Shade B1
- 360034 Shade A2
- 360035 Shade A1
- 360037 Shade C2
- 360038 Shade BL

**ACCESS® TEMPORARY CEMENT**

- 310521 1 x 5 mL Automix Cartridge, Mixing Nozzles, 10-ct.

**NOMIX® TEMPORARY CEMENT**

- 310520 Starter Kit, 1 mL prefilled syringe, 6-ct., Take-Home Applicator, .50 g, 24-ct., 30 Syringe Tips
- 310157 Take-Home Applicator, .50 g, 50-ct.
- 290040 Syringe Refill Kit, 12 x 1 mL prefilled syringe, 60 Syringe Tips
- 640002 16 ga Luer Lock Needle Tip, 30-ct.

DELIVERY SYSTEMS

**ACCUDOSE® & NEEDLETUBE®**

- 290032 19 ga NeedleTube, 100-ct.
- 220008 White Plug, 100-ct.

- 290031 20 ga NeedleTube, 100-ct.
- 220000 Green Plug, 100-ct.

- 290049 24 ga NeedleTube, 100-ct.

- 920201 AccuDose TriPack Starter Kit, 75 Tubes & Plugs

- 210021 Low-Viscosity Tube, 100-ct.
- 210023 Low-Viscosity Tube, 500-ct.
- 210000 Gray Plug, 100-ct.
- 280021 High-Viscosity Tube, 100-ct.
- 280023 High-Viscosity Tube, 500-ct.
- 290046 Black Plug, 100-ct.
- 210001 Separate Clear Tube, 100-ct.
- 210003 Separate Clear Tube, 500-ct.
- 240001 Separate Black Tube, 100-ct.
- 250011 PCR Clear Tube, 100-ct.
- 250001 PCR Opaque Tube, 100-ct.
- 310032 22 ga Luer Lock Needle Tip, 20-ct.
- 640001 22 ga Luer Lock Needle Tip, 100-ct. for Onyx & Champ Syringes
- 640002 16 ga Luer Lock Needle Tip, 30-ct. for Access Edge & NoMIX
- 290038 22 ga Needle Tip, Gray, 144-ct.

**SUPERMIXER™ MIXING NOZZLES**

- 910628 Yellow hub, 50-ct.
- 910631 Yellow hub, 18 ga tip, 50-ct.
- 910632 Pink hub, 50-ct.
- 360010 Teal hub, 24-ct.
- 910626 Blue hub, 50-ct.
- 910630 AutoMix 1:1 Mixing Nozzles, 10-ct.

**BENDA® BRUSH – BRUSH APPLICATORS**

- 370101 Regular, Assorted, 144-ct.
- 370501 Regular, Red, 576-ct.
- 370502 Regular, Blue, 576-ct.
- 370503 Regular, Yellow, 576-ct.
- 370504 Regular, Green, 576-ct.

- 370111 Mini Assorted, 144-ct.
- 370601 Mini, Red, 576-ct.
- 370602 Mini, Blue, 576-ct.
- 370603 Mini, Yellow, 576-ct.
- 370604 Mini, Green, 576-ct.

- 370105 Twin, Pink, 200-ct.
- 370106 Twin, Blue, 200-ct.
- 370107 Twin, Pink, 400-ct.
- 370108 Twin, Blue, 400-ct.

**BENDA® MICRO – FIBER APPLICATORS**

- 380005 Micro Regular, Assorted, 400-ct.
- 380014 Micro Regular, Blue, 576-ct.
- 380015 Micro Regular, Green, 576-ct.
- 380016 Micro Regular, Purple, 576-ct.
- 380017 Micro Regular, Silver, 576-ct.

- 380006 Micro Fine, Assorted, 400-ct.
- 380012 Micro Fine, Yellow, 576-ct.
- 380013 Micro Fine, Red, 576-ct.

- 380033 Micro XF, Assorted, 400-ct.
- 380029 Micro XF, Light Green, 576-ct.
- 380030 Micro XF, Ivory, 576-ct.

- 380021 Micro Twin Regular, Blue, 200-ct.
- 380023 Micro Twin Regular, Blue, 400-ct.
- 380020 Micro Twin Regular, Pink, 200-ct.
- 380022 Micro Twin Regular, Pink, 400-ct.
- 380028 Micro Twin XF, Green, 400-ct.
- 380027 Micro Twin XF, Ivory, 400-ct.



FREE LUNCH & LEARN

- Easy. Esthetic. Durable. Temporization.
- Prevention for Life®
- Impressioning
- The Core Issue

To schedule a Centrix Lunch & Learn, call 800.235.5862 or visit [centrixdentalllearning.com](http://centrixdentalllearning.com)

CENTRIX INCORPORATED

770 River Road | Shelton, CT | 06484 | USA  
203.929.5582, Toll-Free (USA & Canada): 800.235.5862

Centrix is an ISO 13485 certified dental manufacturer and our products are marketed world-wide with FDA, CE, Health Canada, TGA and PAL registrations on file. © 2021 Centrix, Inc. | #GUIDES, REV D

