



Infection Control Offers
for Dental Professionals

2024 PROMOTIONS

June to July 2024

DirectaDentalGroup

DIRECTA topdental **ORSING** parkell **KOHLER** RÖNVIG
TrollDental Sendoline **physics** **CONTACTEZ** ceramir

V-WIPE ZERO

Surface Disinfectant Wipe

BIODEGRADABLE†

Bossklein ECO-XL wipes are constructed from 100% natural plant based fibres. The 40gsm, viscose wipe material is made from 100% renewable resources and approx twice as thick as most competitor wipes.

Each pack contains 100 large alcohol and plastic free wipes each measuring 240 x 200mm.



ECO-XL

100 PLASTIC FREE WIPES
240 x 200mm

EN14476
VIRUS TESTED

2in1
DISINFECTS
& CLEANS

STAY FRESH LID

CONTACT 60 SEC*

APPLE AROMA

PLANT BASED WIPE



Plastic free wipe
Plant based fibres
Renewable resources
Alcohol free
Pack of 100 wipes

BOS8001 V-WIPE ZERO ECO-XL
£8.00 per pack

BUY 9 PACKS FOR £6.00 each BUY 18 PACKS FOR £5.60 each

† Wipes are made from biodegradable plant fibres
*Always test an inconspicuous area first. See website for full details.
Do not use on copper or aluminium.

IMPactiv | P

Impression Disinfectant

Dental impressions are easily contaminated with pathogenic organisms due to their very nature - they come in contact directly with the inside of patient's mouths. Decontamination and disinfection is vital to avoid not only patient/staff cross contamination during application, transport and processing, also to protect the surface quality of the impression itself.



EN14476
VIRUS TESTED

2%
DILUTION
20g/1L

CONTACT 60 SEC

POWDER BASED

DOSING SCOOP

SAFE FOR ALL TYPES

BOS9823 IMPactiv Powder

ONLY £19.95 each

700G

V-WIPE ZERO

Surface Disinfectant Wipe

ALCOHOL FREE

- Refillable dispenser tubs
- Economic refill packs
- Subtle apple aroma
- Effective against Coronavirus
- Safe to use on most surfaces*
- 2 in 1: Contains detergent

* Always test an inconspicuous area first

EN14476
VIRUS TESTED

2in1
CLEANS &
DISINFECTS



CONTACT
60 SEC



APPLE
AROMA



WIPES &
SPRAYS

BOS3000 V-WIPE ZERO Ultra Tub 160 - £11.60 per tub

BUY 6 **FOR £8.70 each** BUY 12 **FOR £8.10 each**

BOS3001 V-WIPE ZERO Ultra Refill 160 - £8.65 per pack

BUY 12 **FOR £6.50 each** BUY 24 **FOR £6.05 each**

BOS2000 V-WIPE ZERO Classic Tub 200 - £9.70 per tub

BUY 12 **FOR £4.95 each**

BOS2001 V-WIPE ZERO Classic Refill 200 - £7.50 per pack

BUY 12 **FOR £5.60 each**

BUY 24 **FOR £5.25 each** BUY 48 **FOR £3.95 each**



40GSM
DURABLE WIPE

18GSM
EVERYDAY WIPE

ULTRA

160 PREMIUM WIPES
175 x 155mm

CLASSIC

200 LARGE WIPES
200 x 200mm

V-WIPE


Surface Disinfectant Wipe

CONTAINS ALCOHOL


- Contact time of only 60 secs
- Wide effective spectrum
- Quick drying
- Leaves no smears or residue
- Dual registration; Medical Device and Biocide

EN14476
VIRUS TESTED


63%
ALCOHOL
CONTENT



CONTACT
60 SEC



LEMON
AROMA



WIPES &
SPRAYS

BOS1000 V-WIPE Ultra Tub 160 - £13.95 per tub

BUY 6 **FOR £10.45 each** BUY 12 **FOR £9.75 each**

BOS1001 V-WIPE Ultra Refill 160 - £11.35 per pack

BUY 12 **FOR £8.55 each** BUY 24 **FOR £7.95 each**

BOS9783 V-WIPE Classic Tub 200 - £12.70 per tub

BUY 12 **FOR £6.25 each**

BOS9784 V-WIPE Classic Refill 200 - £10.40 per pack

BUY 12 **FOR £7.80 each**

BUY 24 **FOR £7.25 each** BUY 48 **FOR £4.95 each**



40GSM
DURABLE WIPE

18GSM
EVERYDAY WIPE

DAILYasp

Aspirator Cleaner & Disinfectant

For the daily cleaning and disinfection of all types of aspirator/suction systems. Can be used on wet line, dry line and semi wet/dry line systems

- Low foaming and non-abrasive
- Low hazard classification
- EN14476 Virus tested
- Concentrate solution
- 2% - 5L makes 250L of solution

EN14476
VIRUS TESTED

2%
DILUTION
20ml/1L

CONTACT
15 MIN

LOW
FOAMING

SAFE FOR
ALL TYPES

100%
CHLORINE
FREE



1 LITRE



2.5 LITRE



5 LITRE

BOS2487 DAILYasp 1 Litre
ONLY **£9.95 each**

BOS2489 DAILYasp 2.5 Litre
ONLY **£24.95 each**

BOS2488 DAILYasp 5 Litre
ONLY **£39.95 each**

UNIoil

Universal Handpiece Oil

All purpose handpiece lubricant spray for maintenance of all dental handpieces.

- Effective lubrication
- Aerosol based
- Peppermint aroma
- Removes debris



BOS9777 UNIoil Handpiece Oil Spray
ONLY **£10.95 each**

easytips

3-in-1 Syringe Tips

SINGLE USE


- Air/water
- Fits UK style syringes
- Use with adapters
- Contains 250 tips

7 CORE
DESIGN



BEVELLED
EDGES

SECURE
FITTING



EAS0001 Easy Tips x 250
ONLY **£19.95 each**