

PROMOTIONS

1ST OF JANUARY - 31ST OF MARCH

PANAVIA™ V5 INTRO KIT

SPECIAL
PRICE

£109.95*

RRP £ 147.66



PANAVIA™ V5, INTRODUCTORY KIT UNIVERSAL (A2) OR CLEAR

1 automix syringe (2.4 ml), 1 PANAVIA™ V5 Tooth Primer (2.0 ml),
1 CLEARFIL™ CERAMIC PRIMER PLUS (2.0 ml), 10 mixing tips,
50 application brushes (fine), 1 mixing dish

*Recommended retail price £ 147.66

**MORE
INFORMATION**



BORN IN JAPAN

Kuraray Europe Benelux, UK & Ireland
Margadantstraat 20, 1976 DN IJmuiden, The Netherlands, info@kireland@kuraray.com, kuraraynoritake.eu
Fiona Bell, mobile: +44 (0)7730 219 605, e-mail: fiona.bell@kuraray.com
Grant Turnbull, mobile: +44 (0)7818 457 568, e-mail: grant.turnbull@kuraray.com

PANAVIA™ V5



ONE CEMENT. ALL CEMENT INDICATIONS. ONE PRIME PROCEDURE.

PANAVIA™ V5 is the cement for all cement indications, irrespective of the indication or the material to be bonded! The recognisable and predictable procedure of PANAVIA™ V5 will take you by surprise with its unprecedented simplicity. Always the same cementation procedure; the special PANAVIA™ V5 Tooth Primer to pre-treat the dental tissue and CLEARFIL™ CERAMIC PRIMER PLUS to prime the restoration. PANAVIA™ V5 is the strongest cement there is! This in combination with no less than five aesthetic colours. The amine-free paste strengthens the durability and colour stability of PANAVIA™ V5.

THE SECRET BEHIND PANAVIA™ V5

The original and purest MDP monomer of Kuraray Noritake Dental still provides the most durable compound with HAp (hydroxyapatite) but also strongly bonds to metals and zirconia. The MDP monomer in PANAVIA™ V5 Tooth Primer and CLEARFIL™ CERAMIC PRIMER PLUS ensures exceptionally durable bonding. This is one of the secrets behind the outstanding performance of PANAVIA™ V5!



INDICATIONS

Cementing of:

- Crowns, bridges, inlays, and onlays
- Veneers
- Adhesion bridges and splints
- Prosthetic restoration on implant abutments and frames
- Posts and cores
- Amalgam bonding

FEATURES

- Dual-cure:
- Automix syringe
- 5 aesthetic colours
- Fluoride release
- Amine-free
- Simple excess removal
- Shelf life of 3 years



Strategic Growth Manager
 Fiona Bell
 Mobile: +44 (0)7730 219 605
 E-mail: fiona.bell@kuraray.com



Business Development Manager
 Grant Turnbull
 Mobile: +44 (0)7818 457 568
 E-mail: grant.turnbull@kuraray.com