

#whdentalwerk



video.wh.com



Show casing W&H

Special offers
January - June 2025



New
Lisa Mini

Speeds you up!

Restoration & Prosthetics

Contra-angles

High torque and precision

1:5 Speed increasing contra-angles, with FG turbine burs, are recommended for cavity and crown preparation. Due to the maximum speed of 200,000 rpm, aerosol creation is reduced.

1:1 contra-angles are recommended for finishing and polishing.



**Synea Fusion
WG-99 LT / WG-99 A**

- › Contra-angle 1:5
- › Optic / Without light
- › 4x Spray
- › Head Ø 9.5 mm

£ **695**

RRP £1,463

£ **628**

RRP £1,308



**Synea Fusion
WG-56 LT / WG-56 A**

- › Contra-angle 1:1
- › Optic / Without light
- › Single Spray
- › Head Ø 9.5 mm

£ **579**

RRP £1,157

£ **462**

RRP £923



Alegria WE-56

- › Contra-angle 1:1
- › Without light
- › Optional External Spray
- › Head Ø 10.0 mm

£ **135**

RRP £356



Alegria HE-43

- › Straight Handpiece 1:1
- › Without light
- › Optional External Spray

£ **125**

RRP £309

Contra-angle Synea Power Edition

For the hardest materials

The Synea Power Edition is especially designed for heavy-duty applications, such as the removal of decayed materials, cavities, fillings, crowns, and bridges. They are the optimal solution for advanced ceramic material like **zirconia** that is increasingly being used in restorative dentistry.



**Launch
Price
55% Off**

**Synea Fusion
WG-900 LT / WG-900 A**

- › Contra-angle 1:4
- › Optic / Without light
- › 3x Spray
- › Head Ø 10.7 mm

£ **695**

RRP £1,489

£ **627**

RRP £1,334



**Synea Vision
WK-900 LT**

- › Contra-angle 1:4
- › Optic
- › 3x Spray
- › Head Ø 10.7 mm

£ **875**

RRP £1,750



Turbines

High speed, low weight

All W&H turbines have a hygienic head* anti-retraction system to prevent cross infection.



**Alegria
TE-97/98**

- › Without light
- › 3x Spray
- › Head Ø 10.4 / Ø 12.2 mm
- › 16 / 18 W
- › 390,000 / 330,000 rpm
- › Roto Quick® Connection

£ **245**

RRP £568



**Synea Vision
TK-94 L**

- › LED+
- › 5x Spray
- › MICRO Head Ø 9 mm
- › 14 W/410,000 rpm
- › Roto Quick® Connection

£ **751**

RRP £1,502



**Synea Fusion
TG-97 / TG-97 L**

- › Without light / LED+
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

£ **503**

RRP £1,006

£ **579**

RRP £1,157



**Synea Fusion
TG-98 / TG-98 L**

- › Without light / LED+
- › 4x Spray
- › MIDI Head Ø 11.5 mm
- › 21 W/360,000 rpm
- › Roto Quick® Connection

**Special Offer Price
When you buy
3x TG 97/98**

**£ 1,358 Non-optic
£ 1,562 Optic**

Prophylaxis

Polishing by rotation

Enjoy professional tooth cleaning and give your patients a beaming smile.



**Proxeo
HP-44 M**

- › Handpiece 4:1
- › System Doriot
- › For disposable contra-angles with plastic shank

£ **199**

RRP £328



**Proxeo
WP-64 MU**

- › Prophy contra-angle 4:1
- › Without light
- › For universal screw-in and snap on (adaptor) cups and brushes

£ **289**

RRP £578

Air motors & Couplings



Couplings
Available from

£ **105**

RRP £210



Air motors
Available from

£ **247**

RRP £448



Our turbines are available to fit other manufacturers couplings.

Oral Surgery & Implantology

Implantmed Plus SI-1023

Workflow driven solution

Flexibility to meet your needs!

W&H offer a synergistic solution for oral surgery guiding the clinician through the cycle from implementation to evaluation.

You can enjoy the benefits of utilising the modular Implantmed surgical motor and Piezomed piezosurgery unit alongside the Osstell Connect consultation tool and the Osstell Beacon, which ensures predictable assessment of ideal loading.




Free Surgical WS-75 Contra-angle Included

with selected Implantmeds

Piezomed Module

Plug & Play piezo surgery device for Implantmed Plus

Display menu and documentation in your local language



Motors EM-19 EM-19 LC

With and without light

Osstell Beacon

Measure and evaluate the stability of the implant

Straight & Contra-angle Handpieces

Wireless foot control for controlling multiple W&H devices



Implantmed Plus SI-1023 Kit 01

Non-optic/Cable

- Non-optic motor
- Cable Foot Control

FREE WS-75

£ 3,690

RRP £7,379



Implantmed Plus SI-1023 Kit 04

Optic/Wireless

- Optic motor
- Wireless Foot Control

FREE WS-75 L

£ 4,610

RRP £9,220



LIGHT



NO LIGHT



CABLE FOOT CONTROL



WIRELESS FOOT CONTROL

2 Years* Warranty

Implantmed Classic SI-923

Implantmed excels by virtue of its ease of operation, powerful motor and automatic thread cutter function.



Implantmed Classic SI-923

Non-optic/Cable

- Non-optic motor
- Cable Foot Control

FREE WI-75 E/KM

£ 3,200

RRP £6,401

Oral Surgery & Implantology

Piezomed

"A real practice builder"



Piezomed Classic Optic/Cable

- Piezomed Classic
- Handpiece with LED ring
- Cable Foot Control
- FREE BONE SET Worth £719

£ 5,310

RRP £9,671



Piezomed Plus Module

- Power 24W. Handpiece with light SA-40 L
- Automatic instrument detection
- Instrument changer and instrument library
- FREE EXTRACTION SET INCLUDED

£ 3,411

RRP £6,201

Piezomed Instruments

The W&H Piezomed instrument range includes different instruments for bone surgery, periodontology and extraction, for retrograde endodontics and for lateral sinus lift.



Automatic instrument detection



Bone surgery

£ 741

RRP £1,348



Extraction

£ 235

RRP £428



Periodontology

£ 328

RRP £596



Retrograde endodontics

£ 582

RRP £1,058



Lateral sinus lift

£ 532

RRP £967



Preparation of the implant site and crestal sinus lift

£ 1,043

RRP £1,897

Oral Surgery & Implantology

Surgical Handpieces

High-quality stainless steel with a **scratch-resistant coating** makes W&H surgical handpieces robust and gives them a particularly long service life.



Self-generating LED Surgical Handpieces available.



Contra-angle WS-91 / 1:2.7

- › Without light
- › 45° head with 3x spray
- › Dismantleable

£ **757**

RRP £1,513



Contra-angle 20:1 WI-75 E/KM

- › Without light
- › External single spray
- › Not dismantleable

£ **650**

RRP £1,300



Contra-angle WS-91 L / 1:2.7

- › Mini LED+
- › 45° head with 3x spray
- › Dismantleable

£ **837**

RRP £1,674



Contra-angle 20:1 WS-75

- › Without light
- › External single spray
- › Dismantleable

£ **803**

RRP £1,606



Straight handpiece S-11 / 1:1

- › Without light
- › External single spray
- › Dismantleable

£ **649**

RRP £1,297



Contra-angle 20:1 WS-75 L

- › Mini LED+
- › External single spray
- › Dismantleable

£ **879**

RRP £1,757



Straight handpiece S-11 L / 1:1

- › Mini LED+
- › External single spray
- › Dismantleable

£ **913**

RRP £1,826



Handpiece 20:1 Zygoma SZ-75

- › Without light
- › External single spray
- › Dismantleable

£ **1,190**

RRP £2,527



Oral Surgery & Implantology

Osstell Beacon

Measurement of implant stability




Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the easy to read colour light indicator and ISQ reading. Reduce treatment times and avoid rejection in high-risk patients.

£ **1,995**

RRP £5,277




Osstell SmartPegs

The SmartPeg measures the resonance frequency when it is attached to an implant.

£ **85**

RRP £117



What are the positive outcomes when using the Osstell Implant stability solutions?

- › A positive experience for patients, encouraging trust through understanding
- › Enhanced patient education builds confidence, resulting in increased treatment acceptance
- › Reduction in treatment visits and failures with improved average healing times
- › Improved outcomes for patients with risk factors
- › Enhanced communication between specialists ensuring a more synergistic approach to treatment
- › Improved treatment documentation

Osstell ISQ devices with the Osstell ISQ Scale help implantologists to make informed implant treatment decisions.



Osstell Connect FREE App for Osstell Beacon

Osstell Connect is designed to document quantified data which ensures traceability, whilst complimenting data in your Patient Management System.

Systems with graphical presentation of data and statistics which is great for patient communication. You can easily present treatments in a graphic and totally personalized way, even generating pdf reports, while having access to the world's largest online surgical portal for peace of mind. The patient consultation tool allows you to easily communicate your proposed implant.

*Register online at Osstell Connect to claim your free SmartPegs when you buy an Osstell Beacon



New!

BePro Cleaner TE

Tri-enzymatic non-foaming detergent for automatic cleaning in washer disinfectors.

£ 45

RRP £56



BePro Rinse

For rapid, stain-free drying during automated reprocessing of dental instruments.

£ 14

RRP £18

Teon thermal washer disinfectors

- › Validated cleaning of all instruments
- › Traceability and connectivity through Ethernet module
- › 1 Year Warranty

Teon with Standard Starter Set

£ 6,875

RRP £9,365

Teon with Handpiece Injector Starter Set

£ 8,200

RRP £10,590

BePro WDD Starter Pack

£ 72

- › BePro WDD Cleaner TE – 5000 ml
- › BePro WDD Rinse – 1000 ml
- › Canister opener



134°C
SSS

Sterilization at its Best!

With a wide range of sterilizers, W&H is able to offer the perfect option to meet individual needs.



S 23'

Lyla S type Sterilizer New concept sterilizer

- › Mobile App available via upgrade for remote data storage
- › Future-proof upgradability to a type B sterilizer using the new Activation Code concept

Unique upgrade flexibility with **Activation Codes**

- › B Cycle
- › Remote data storage
- › Traceability
- › All in One

£ 3,400

17 l. RRP £4,762

£ 3,800

22 l. RRP £5,008

*Based on average 2kg load for immediate use



B 20'

Lara B type Sterilizer

- › Mobile App for Real Time Remote Monitoring
- › Eco Dry technology automatically adapts the drying time to the mass of the load
- › Unique optional upgradability to increase: Performance and Traceability or to add a Fast Cycle

Unique upgrade flexibility with **Activation Codes**

- › Performance**
- › Traceability
- › Fast Cycle
- › All in One

** Only for Lara 17 and 22 litres

£ 4,650

17 l. RRP £6,857

£ 5,250

22 l. RRP £7,700

*with fast cycle upgrade

Lara XL

£ 5,750

28 l. RRP £8,492



B 13'

Lisa B type sterilizer

- › Full extended traceability thanks to EliTrace
- › Mobile App for Real Time Remote Monitoring
- › Eco Dry+ technology automatically adapts the drying time to the mass of the load
- › Elisense user interface generates status, temperature and other cycle information
- › Cycle times from 13 minutes for specific loads
- › State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use

£ 5,550

17 l. RRP £8,467

£ 6,290

22 l. RRP £9,580



B 10'

NEW Lisa MINI B type sterilizer

- › Cycle times from 10 minutes
- › Integrated vacuum pump
- › Small and easily stackable: W 44cm x H 20cm x D 67cm
- › Built in water tank and automatic water filling - no external storage
- › For up to 22 instruments, offers complete traceability
- › 5 l. model includes internal vacuum pump
- › Includes integrated water tank

3 liters version available
Requires compressed air

£ 3,400

5 l. RRP £4,900



Sterilization, Hygiene & Maintenance

All W&H Sterilizers offer:

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Cycle traceability available thanks to Wi-fi connection and built in data logger
- › Ergonomic design, with modular feet for flexible positioning
- › Easy to clean surfaces and accessible water › Automatic Water Fill function
- › Pre-programmable vacuum and helix (for B Type Sterilizers) test function to run automatically out of surgery hours

Sterilization, Hygiene & Maintenance

Internal cleaning & lubrication

Increase the life of your handpieces and reduce the costs of repairs!



10" quick connect by W&H

Assistina Twin

- > Two process chambers allowing loading of a 2nd handpiece while the 1st is processed
- > Superior Maintenance & Lubrication
- > Intuitive and easy to use
- > A flexible range of adaptors available

£ 1,441

RRP £3,202



18" quick connect by W&H

Assistina One

- Internal cleaning & lubrication of one handpiece in only 18 seconds, with Assistina One the aerosols stay in the HEPA filter thereby releasing clean air.

£ 1,020

RRP £2,169

Assistina Twin and Assistina One includes

Care set and Quick RM and ISO adaptors



Service Oil F1

Ask your dealer about W&H Handpiece Oil



3x EliReel Included

Seal2 Starter Pack

- > Universal Pouch Sealing
- > 12mm hermetic seal
- > 3 x EliReel Rolls FREE

£ 747

RRP £1,245

W&H EliReel sealing pouches

From / per roll

£ 14

RRP £22



Osmo

Supplies multiple sterilizers simultaneously with top quality demineralized water.

- > Automatic water filling
- > Most cost-effective solution
- > Increase sterilizer service life
- > Up to 6 devices simultaneously

£ 926

RRP £1,543



Seal Check

Pack of 100 sealing check

£ 31

RRP £62



Lisa Safe Label Printer & QR Barcode Reader

Quick and easy tool for linking the sterilization cycle and dental products used for treatment

Barcode reader

£ 145

RRP £234

Lisa Safe printer

£ 510

RRP £802

Pro Service

PremiumCare



Gives you the flexibility to select a plan that best suits the individual needs of your practice.

The standard PremiumCare package offers a flexible plan allowing you to add optional extras as you please or the **PremiumCare Plus or Flex**** plan, which include fully comprehensive service packages and the required pressure vessel inspection and validation testing at no extra cost.

PremiumCare

Maintenance plan for all W&H sterilizer models, with a maximum limit of 4,000 cycles over the 3 year term*.

PremiumCare Plus & Flex**

A fully comprehensive maintenance plan for all W&H sterilizer models, including PVI certification, unlimited cycles and validation in line with manufacturer's or current HTM 01-05 guidance (model dependent).

Our Technical Service Team

are here to support you nationwide!




ProService

Services	Premium Care	Premium Care Plus & Flex**
All routine service visits for three years	✓	✓
All service parts and scheduled sundries	✓	✓
All breakdown calls	✓	✓
All breakdown parts including labour	✓	✓
All engineer travel in the UK Excluding Highlands, Islands and ferries etc,	✓	✓
Free on-site operator training	✓	✓
Free technical software upgrades	✓	✓
Unlimited Cycles - for the duration of the plan	⊖	✓
PVI Certification - in line with Pressure Systems Safety Regulations 2000:128	⊖	✓
Periodic Test Logbook	⊖	✓
Validation in line with manufacturer's or current HTM 01-05 guidance (model dependent)	⊗	✓
Installation and Commissioning in line with current guidelines	⊗	✓
Helix Test Kit***	⊗	✓
Extra Discount	⊗	✓



Included



Not Included



Optional

BePro

Professional Cleaning and
Disinfection for All Surfaces

be pro



Soaking bath

Soaking bath for the immersion
of dental instruments in the
decontamination phase.

£ 38

RRP £46



Disinfectant R 1L.

Ready-to-use solution for the
disinfection and cleaning of
rotary dental instruments.

£ 13

RRP £16



Disinfectant I 2L.

Powerful disinfectant
concentrate for medical
and dental instruments.

£ 38

RRP £48



Disinfectant Spray 1L.

Surface disinfectant for the
rapid disinfection and cleaning
of non-invasive medical devices.

£ 13

RRP £16



Disinfectant Wipes

120 wipes for the
disinfection of non-
invasive medical devices.

£ 9

Tube RRP £11

£ 7

Refill RRP £9

New!

Meet the Sales Team



Jon Bryant

National Sales and Marketing Manager

07917 886326

jon.bryant@wh.com



Richard Mansbridge

Area: London, South East, Anglia,
East Midlands and Humberside

07775 792628

richard.mansbridge@wh.com



Caroline Caine

Area: West Midlands, Manchester, Merseyside
South West England, Isle of Man, Channel
Islands, Wales, Northern Ireland and Eire

07786 515331

caroline.caine@wh.com



David Barrass

Area: Scotland, Humberside,
North of England and North West

07795 315089

david.barrass@wh.com

W&H (UK) Ltd

6 Stroud Wood Business Centre,
Park Street, St Albans, AL2 2NJ

t +44 (0) 1727 874990

f +44 (0) 1727 874628

office.uk@wh.com

www.wh.com

To place an order
contact W&H UK on
+44 (0) 1727 874990
or contact your preferred
dental dealer

Offers valid from 1st January to 30th June 2025. Information is correct at time of publication and relates to W&H products and services available in the UK and Ireland. Prices are quoted in £ sterling and exclude VAT.

Savings are based on RRP and rounded to the nearest £.

Free goods are not transferrable and no cash equivalent is available.

Terms and Conditions apply. E. & O. E.