




ortho

r  
t  
h  
o  
d  
o  
n  
t  
i  
c  
s



**Trycare** for all your **Aligner Solutions**



### Dentalastics Separators




Blue for Posteriors. For quick and easy separation of teeth. Radiopaque. Inner diameter 2.1mm.

Blue	Pack 1000	UWDS000
------	-----------	---------



### Elastic Separators With Stops



Safe-Seeps Elastomeric Separators, Regular Blue.

Blue	Pack 1000	UPRO009
------	-----------	---------



### Separating Plier



Ideal for the placement of elastic separators.

Each	UHFPOSP
------	---------




### IPR Gauge



Includes 6 Stainless Steel gauges of 0.1 to 0.5mm.

UFSIPRG
---------




Dr. Francesco Garino, Turin, Italy



### IPR Distance Control Set

- Stainless Steel and autoclavable.
- A complete and practical IPR Gauge.
- Safe grip thumb & finger holder.
- 8 sizes correlate to the Intensiv IPR range.  
8 Sizes: 0.10mm to 1.0mm.

UITV006
---------



### ApproxOpener

Designed for opening of contact points with a serrated strip edge and 8µm diamond coating for initial IPR. Safe grip thumb & finger holder.

		Thickness	
Double-Sided	Pack 3	0.09mm	UITV001
Single-Sided	Pack 3	0.10mm	UITV002



### Lightning Strips Medium

Size	Single Sided	Pack	Code	Size	Double Sided	Pack	Code
4mm	Single Sided	12	UGHF457	4mm	Double Sided	12	UGHF460
6mm	Single Sided	12	UGHF458	6mm	Double Sided	12	UGHF461
8mm	Single Sided	12	UGHF459	8mm	Double Sided	12	UGHF462



### Separating Strips

Width	Single Sided	Thickness	Code
3mm	F303	0.07mm	UPRO004
4mm	C304	0.09mm	UPRO005
6mm	C306	0.09mm	UPRO006
6mm	416 Serrated	0.05mm	UPRO007




Corundum-bound, electro-deposited Stahlcarbo® strips. Can be cut to desired length. Packet of 12 strips.



### Diamond Finishing Strips - Perforated - Single Sided

Gain space by gently contouring the sides of the overlapping teeth. Made of high grade Stainless Steel to resist breaking and stretching. Autoclavable. Available in a variety of widths and grain.



#### Perforated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0.08mm	UEDFS18
Fine	3.75mm	30µm	0.10mm	UEDFS17
Medium	3.75mm	45µm	0.13mm	UEDFS16


#### Perforated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0.08mm	UEDFS12
Fine	2.5mm	30µm	0.10mm	UEDFS13
Medium	2.5mm	45µm	0.13mm	UEDFS11



### Diamond Finishing Strips - Perforated Serrated - Single Sided

A serrated saw edge and diamond coating combined in one strip, allowing for simpler working procedures.



#### Perforated Pack 10


Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0.08mm	UEDFS09
Fine	3.75mm	30µm	0.10mm	UEDFS08
Medium	3.75mm	45µm	0.13mm	UEDFS07

#### Perforated Serrated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0.08mm	UEDFS15
Fine	2.5mm	30µm	0.10mm	UEDFS14
Medium	2.5mm	45µm	0.13mm	UEDFS10



### Diamond Finishing Strips - Plain Serrated - Single Sided



#### Perforated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0.08mm	EDNT0265


#### Perforated Serrated Pack 10

Colour	Width	Grain	Thickness	Code
Fine	3.75mm	30µm	0.10mm	UEDFS04



### Diamond Finishing Strips - Plain - Single Sided

Diamond Interproximal Strips ideal for enamel reduction as recommended for aligner treatment. High grade Stainless Steel to resist breaking and stretching. Autoclavable.



Pack 10					Pack 10				
Colour	Width	Grain	Thickness	Code	Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0.08mm	UEDFS03	Extra Fine	2.5mm	15µm	0.08mm	UEDFS02
Fine	3.75mm	30µm	0.10mm	UEDFS05	Fine	2.5mm	30µm	0.10mm	UEDFS01
Medium	3.75mm	45µm	0.13mm	UEDFS06	Medium	2.5mm	45µm	0.13mm	UEDFS00



### Diamond Finishing Strips - Plain Wide - Single Sided



Pack 5				
Colour	Width	Grain	Thickness	Code
Extra Fine	6mm	15µm	0.08mm	UEDCDS6
Fine	6mm	30µm	0.10mm	UEDFDS6
Medium	6mm	45µm	0.13mm	UEDFS19

#### Essix

### Flex-View™ Single Sided

These perforated Stainless Steel strips allow optimal visibility and minimise clogging. Flex-View™ Interproximal Strips also feature a gateway area for easy access for interproximal stripping and finishing. All strips are single sided and autoclavable.



#### 2.5 mm Narrow Strips - Pack 10

X-Fine Grit Yellow	UORC095
Fine Grit Red	UORC093
Med Grit Blue	UORC094

#### Pack 15

5 x Wide Extra Fine, 5 x Fine, 5 x Medium	Assorted	UORC092
---	----------	---------

#### 3.75 mm Wide Strips - Pack 10

X-Fine Grit Yellow	UORC099
Fine Grit Red	UORC097
Med Grit Blue	UORC098


#### Pack 15

5 x Wide Extra Fine, 5 x Fine, 5 x Medium	Assorted	UORC096
---	----------	---------

**Essix**

**Diamond Interproximal Strips**

Fully autoclavable, surgical steel plain strips. These very thin strips (0.08mm/0.003") are made of electrolytically bonded diamonds on a vibro cushion base.



**Single Sided Pack 10**

Colour	Width	Grain	Thickness	Code
Fine	2.5mm	30µm	0.08mm	UORC107
Fine	3.75mm	30µm	0.08mm	UORC110
Medium	2.5mm	45µm	0.08mm	UORC106
Medium	3.75mm	45µm	0.08mm	UORC109

**Double Sided Pack 10**

Colour	Width	Grain	Thickness	Code
Medium	2.5mm	45µm	0.08mm	UORC111
Medium	3.75mm	45µm	0.08mm	UORC108




**Essix**

**Interproximal Serrated Strips 6mm**

Stainless steel ultra thin interproximal serrated strips, used prior to reduction with diamond strips. Width: 6mm, Thickness: 0.05mm.

Pack 10	UORC105
---------	---------



**Komet**

**Serrated Saw Diamond Strips**

Diamond Strips. 0.05mm Serrated Separating Saw. Width: 6mm, Thickness: 0.05mm.

Pack 10	TD0703
---------	--------



**Komet**

**Vision Flex Diamond Strips**

These Strips have a more dead-soft quality than the typical springy feel of conventional diamond strips.

**Pack 10 - Length 147mm**

Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0.08mm	UVFW37E
Fine	3.75mm	30µm	0.10mm	UVFW37F
Medium	3.75mm	45µm	0.13mm	UVFWS37

**Assorted Pack 15**

Contains: 5 x medium (45µm), 5 x fine (30µm) and 5 x extra-fine (15µm)

Assorted 3.75mm	UVFW37A
-----------------	---------

**Pack 10 - Length 147mm**

Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0.08mm	UVFW25E
Fine	2.5mm	30µm	0.10mm	UVFW25F
Medium	2.5mm	45µm	0.13mm	UVFWS25

**Assorted Pack 15**

Contains: 5 x medium (45µm), 5 x fine (30µm) and 5 x extra-fine (15µm)

Assorted 2.5mm	UVFW25A
----------------	---------



Komet

Solid Diamond Strips



Pack 10 - Length 148mm

Colour	Width	Grain	Thickness	Code
Extra Fine	3.7mm	15µm	0.08mm	UVFD37E
Fine	3.7mm	30µm	0.10mm	UVFD37F
Medium	3.7mm	45µm	0.13mm	UVFDS37
Medium	3.7mm	45µm	0.21mm	UVFD37D
Coarse	3.7mm	90µm	0.18mm	UVFD37C

Pack 10 - Length 148mm

Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0.08mm	UVFD25E
Fine	2.5mm	30µm	0.10mm	UVFD25F
Medium	2.5mm	45µm	0.13mm	UVFDS25



Code: UEDS04K  
PDL-4 KIT 4mm



Code: UEDS25K  
PDL-2.5 KIT 2.5mm



EVE Diastrips

Diastrips (Porcelain, Composite, Metal, Enamel).

- For manual surface processing
- Interdental application
- Perforation removes the cut material.



1 Coarse 2 Medium 3 Fine

PACK 9

Kits Contain: 3 x of each grit (fine, medium and coarse) in the selected size of 2.5mm or 4mm widths.

						<b>Dimensions (mm)</b>
4 x 150	2.5 x 150	4 x 150	2.5 x 150	4 x 150	2.5 x 150	
0.11	0.11	0.10	0.10	0.07	0.07	<b>Thickness (mm)</b>
UEDS454	UEDS452	UEDS354	UEDS352	UEDS104	UEDS102	<b>Code:</b>



## Garrison FitStrip

**FitStrips are your complete line of interproximal trimmers and finishers!**

Whether you're looking for safe and accurate IPR for your aligner cases, interproximal cement clean-up, or smoothing out cervical margins, FitStrip™ is a fast and convenient solution.

FitStrip™ is indicated for composite restorations (contouring/finishing all interproximal surfaces), orthodontics (IPR/tooth slenderizing), and crown and bridge (cement removal and clean-up).



### Garrison FitStrip Starter Kit (FPSK01)

**Starter Kit Includes:**

- 2 x Serrated Strip (0.05mm)

**Single Sided**

- 1 x Super Fine Grit Strip (0.08mm)
- 1 x Fine Grit Strip (0.10mm)
- 1 x Medium Grit Strip (0.13mm)
- 1 x Coarse Grit Strip (0.18mm)

**Double Sided**

- 1 x Super Fine Grit Strip (0.11mm)
- 1 x Fine Grit Strip (0.15mm)
- 1 x Medium Grit Strip (0.21mm)
- 1 x Coarse Grit Strip (0.30mm)
- 2 x FitStrip Handle



Starter Kit UFPSK01



### Garrison FitStrip Universal Kit (FPSK04)

**Universal Kit Includes:**

- 1 x IPR gauge & 2 x FitStrip™ Handle.

**Single Sided**

- 1 x Super Fine Grit Strip (0.08mm)
- 1 x Fine Grit Strip (0.10mm)
- 1 x Medium Grit Strip (0.13mm)
- 1 x Coarse Grit Strip (0.18mm)
- 1 x Super Fine Grit Strip (0.11mm)

**Double Sided**

- 1 x Fine Grit Strip (0.15mm)
- 1 x Medium Grit Strip (0.21mm)
- 1 x Coarse Grit Strip (0.30mm)
- 1 x Coarse Grit Strip (0.40mm)
- 1 x X-Coarse Grit Strip (0.50mm)

Universal Kit UFSUNKT



### Garrison FitStrip IPR Kit (FPSK05)

**IPR Kit Includes:**

- 1 x IPR Gauge
- 2 x FitStrip Handle

**Double Sided**


- 2 x Super Fine Grit Strip (0.11 mm)
- 2 x Medium Grit Strip (0.21mm)
- 2 x Coarse Grit Strip (0.30mm)
- 2 x Coarse Grit Strip (0.40mm)
- 2 x X-Coarse Grit Strip (0.50mm)

Universal Kit UFSIPRK





# Interproximal Finishing & Contouring Strip Refills



## Garrison FitStrip Serrated

FitStrips' simple, detachable handle makes a world of difference for both clinician ease of use and patient comfort.

Your hand stays outside the mouth, granting a clear field of view and easier operation. Twist the colour-coded barrel to adjust FitStrips' curvature and it's automatically locked in place. It doesn't get any easier than that. FitStrips are very durable and long lasting.


Serrated Grey 0.05mm

Pack of 4

UFPSSAW



## Garrison FitStrip



### Double Sided

### Pack Size

### Code:

Superfine Yellow 0.11mm	Pack of 4	UFDSF11
Fine Red 0.15mm	Pack of 4	UFDF015
Medium Blue 0.21mm	Pack of 4	UFPSMDS
Coarse Green 0.30mm	Pack of 4	UFDC030
Coarse White 0.40mm	Pack of 4	UFDMG04
X-Coarse Black 0.50mm	Pack of 4	UFDMB24

### Single Sided

### Pack Size

### Code:

Superfine Yellow 0.08mm	Pack of 4	UFPSFS4
Fine Red 0.10mm	Pack of 4	UFPFSR4
Medium Blue 0.13mm	Pack of 4	UFPSMSS
Coarse Green 0.18mm	Pack of 4	UFPSCS4





### QwikStrip Ultimate Assorted Kit

Industry Leading Qwik Strips.

- Easy to follow IPR sequence
- Interproximal cement removal
- Use for any aligner regime
- Double or Single Sided
- Autoclavable



Scan the QR Code: to watch a QwikStrip IPR Demonstration.

Kit Contains:

- 5 x Single-sided (.07mm - .18mm)
- 5 x Double-sided (.18mm - .40mm)

UQWKS01

Ideal for breaking contacts and for interproximal cement removal.



### QwikStrip Single-Sided

.07mm Yellow	Pack 10	UQWKS03
.10mm Red	Pack 10	UQWKS04
.13mm Blue	Pack 10	UQWKS05
.16mm Green	Pack 10	UQWKS06
.18mm Black	Pack 10	UQWKS07



### QwikStrip Double-Sided

.18mm Yellow	Pack 10	UQWKS13
.20mm Red	Pack 10	UQWKS14
.25mm Blue	Pack 10	UQWKS15
.30mm Green	Pack 10	UQWKS16
.40mm Black	Pack 10	UQWKS17



### QwikStrip Curved


.07mm Yellow	Pack 10	UQWKS08
.10mm Red	Pack 10	UQWKS09
.13mm Blue	Pack 10	UQWKS10
.16mm Green	Pack 10	UQWKS11
.18mm Black	Pack 10	UQWKS12



### QwikStrip Serrated

.05mm White	Pack 10	UQWKS02
-------------	---------	---------





### StripFix Finishing Strip System

StripFIX finishing strip system Manual Diamond Strip Holders. Use with or without the handle. Hold between middle finger with your index finger on the diamond strip plate. Strip holder is swivel mounted. May be locked in 4 positions for optimal access to the oral cavity. Simple quick-plug system to hold the strip. All components autoclavable.



**Pack contains:**

- 1 x .04mm Serrated Edge
- 1 x .10mm Fine
- 1 x .08mm Super Fine
- 1 x .13mm medium

Pack of 4 Files and Handle	USFX000
----------------------------	---------

#### Refills - Single-Sided

White Serrated Edge .04mm	<b>Pack 10</b>	USFX001
Yellow Super Fine .08mm	<b>Pack 10</b>	USFX002
Red Fine .10mm	<b>Pack 10</b>	USFX003
Blue Medium .13mm	<b>Pack 10</b>	USFX004
2 x Handles .13mm	<b>Pack 2</b>	USFX005





**Intensiv**  
Swiss Dental Products  
**Swingle Handpiece**

W&H Synea 2:1 Reciprocating Handpiece for the Intensiv Ortho Strip System. Available with or without light.

With Light	UITV049
Without Light	UITV048

**Clinical Case**



Opening of interdental space



Creating additional space by interproximal enamel reduction (IPR)



Finishing/Polishing.



**Intensiv**  
Swiss Dental Products  
**IPR Set Professional**

Calibrated and precise stripping (IPR) for orthodontics. Mechanical and manual strips in a stainless steel tray. The professional kit contains both manual and mechanical IPR options all in an autoclavable cassette.

REF WG-69A * 1 Piece	Intensiv Swingle 2:1, reciprocating contra-angle (without light)
REF 053 1 Piece	Intensiv Ejector for safe file removal between IPR stages
REF OS08OP-DS 1 Piece	Intensiv Ortho-Strips, 8µm opener, serrated strip edge, for opening of contact points
REF OS40M-DS 1 Piece	Intensiv Ortho-Strips, 40µm medium, for step one creating additional space
REF OS60C-DS 1 Piece	Intensiv Ortho-Strips, 60µm coarse, for step two in creating additional space
REF OS80XC-DS 1 Piece	Intensiv Ortho-Strips, 80µm extra coarse, perforated, for high enamel reduction
REF OS40M-DS 1 Piece	Intensiv Ortho-Strips, 40µm medium, for pre-polishing of the enamel surface
REF OS15POL-DS 1 Piece	Intensiv Ortho-Strips, 15µm polishing, for fine polishing of treated surfaces

All in a single autoclavable box. UITV005

**New  
Central  
System**



**Intensiv OST400Set01/6**

Thanks to the diamond coating in the central operating area, the Intensiv Ortho-Strips System Central allows for high precision in approximal reduction, finishing and polishing of enamel, respecting the tooth morphology and dental tissue.

During the process of enamel reduction, trimming and polishing in orthodontics, parts of enamel or dentine below the tooth margin may be removed unintentionally, which may cause subsequent hypersensitivity of the tooth in the patient.

Colour Code:	Micron	Description
	494	8µm Opener
	504	15µm Polishing
	514	25µm Fine
	514	40µm Medium
	524	60µm Coarse

- Contains:**
- 1 x Ortho-Strips Tray
  - 1 x Ortho-Strips System Opener.

- Ortho-Strips system central double sided:**
- 2 x Contouring
  - 1 x Reduction
  - 1 x Finishing
  - 1 x Polishing.

Kit	UITV010				
Double-Sided diamond					
	Opening	Contouring	Reduction	Finishing	Polishing



**Ortho-Strip Central**

Pack of 3.

Colour Code:	Micron	Description
	524	60µm Coarse
	514	40µm Medium
	514	25µm Fine
	504	15µm Polishing


Double-Sided diamond				
	Reduction	Contouring	Finishing	Polishing
	<b>Code:</b> UITV025	UITV022	UITV019	UITV016







Single-Sided diamond				
	Reduction	Contouring	Finishing	Polishing
	<b>Code:</b> UITV026	UITV023	UITV020	UITV017
	<b>Code:</b> UITV027	UITV024	UITV021	UITV018

**Ejector for Swingle**

Intensiv File Ejector

UITV004









Colour Code:	Micron	Description
	494 8µm	Opener
	504 15µm	Polishing
	514 25µm	Fine
	514 40µm	Medium
	524 60µm	Coarse
	524 80µm	Extra coarse



### Ortho Strips







Fast, efficient and precise, double and single sided.  
Metal frames for reduction and polishing. Pack 3

**Double-Sided diamond**

Code:						
	UITV028	UITV033	UITV034	UITV037	UITV040	UITV045
	Opening	Polishing	Finishing	Contouring	Reduction	Space creation







  

**One-Sided diamond**

Code:						
	UITV029	UITV030	UITV031	UITV032	UITV035	UITV036
	Opening	Opening	Polishing	Polishing	Finishing	Finishing

**One-Sided diamond**

Code:						
	UITV038	UITV039	UITV041	UITV042	UITV043	UITV044
	Contouring	Contouring	Reduction	Reduction	Space creation	Space creation



### Ortho Strips

Contains an assortment of 6 oscillating files:

- Opening 08µm
- Space Creating 80µm
- Polishing 15µm
- Reduction 60µm
- Finishing 25µm
- Contouring 40µm

Double-Sided Set	Each	UITV007
Single-Sided Set - Left	Each	UITV008
Single-Sided Set - Right	Each	UITV009






**Pop On Diamond Finishing Discs  
Double Sided**

Double sided 0.10mm thick.  
Diamond Coated. Each.

10mm	Each	UKMT001
14mm	Each	UKMT002
18mm	Each	UKMT003
Mandrel	Each	UKMT004



**Superflex  
Diamond Open Mesh Disc**

Superflex Diamond Open Meshed Disc  
Double Sided 22mm Diameter, 0.15mm Thickness.

22mm Diameter	Each	USFDO22
---------------	------	---------

**Komet DISCtance IPR system kit**



UKMT005

Rotary diamond discs for interproximal enamel reduction.

**Comes with 5 Disc Types:**

Code	Grit	Thickness	Coating
UWES001	Fine	0.15mm	Upper side
UWES002	Fine	0.15mm	Bottom side
UWES003	Fine	0.15mm	Both sides
UWES004	Medium	0.20mm	Both sides
UWES005	Medium	0.30mm	Both sides

The Komet DISCtance IPR system features a range of new, sharp stripping discs that ensure the amount of enamel is easily, efficiently and precisely removed. All discs are mounted on a straight handpiece shank. Fully autoclavable bur holder.



**IPR Diamond Bur**

IPR Diamond Bur. Pack 5.



**392 Series**

Super Fine	UEI392S
Fine	UEI392F
Medium	UEI392M



**Orthodontic IPR Kit**





**The Komet IPR & Disc kit contains:**

- Vision Flex Strips
- 3x WS37 .13mm 45µm Blue Medium
- 3x WS37F .10 mm 30µm Red Fine
- 3x WS37EF .08 mm 15µm Yellow Extra Fine
- 3x DS37D .21mm 45µm Blue Medium Double Sided
- 3x DS37C .18mm 90µm Green Coarse Single Sided
- 1x 934.104.180 Diamond Disc .18mm 18mm Diameter
- 1x 8934A.900.180 .15mm 14mm Diameter
- 1x Contra Angle Mandrel
- 1x De-Bonding Bur
- 1x H379AGK Egg De-Bonding Bur
- Non-autoclavable storage box

UKMT006

# Precise Adjustments For Better Results



**Hu-Friedy**

## Tear Drop

Creates a notch at the gingival margin of clear aligners anywhere elastic hooks are required within the arch. Provides a reservoir to hold the elastics in place as the aligner is seated, permitting easier manipulation of elastics by the patient.

UHFAPTD

Ref: 678-800



**Hu-Friedy**

## Vertical

Emphasizes over-correction of rotations by producing an an indentation on the facial or lingual aspect of a specific tooth.

UHFAPVT

Ref: 678-801



**Hu-Friedy**


## Horizontal

Enhances retention of either clear aligners or retainers by placing an indent into a tooth's undercut or below a bonded attachment.

UHFAPHZ

Ref: 678-802





**Hu-Friedy**

### Hole Punch

Creates half-moon cut-outs for relief of clear aligners around bonded buttons used for elastics. Nips away plastic to provide soft tissue clearance to prevent aligner impingement of tissue.

UHFAPPN

Ref: 678-803



**Hu-Friedy**

### Petite Punch

The Petite Punch cuts a 4.5mm half-moon cut-out in an aligner for relief of the plastic around bonded buttons used for elastics.

UHFAPP

Ref: 678-805




**Hu-Friedy**

### Spot

The Spot creates a tiny rounded dimple for spot contacts to accent small individual tooth movements or added retention.

UHFAPSP

Ref: 678-804



**Hu-Friedy**

### Counter Clockwise Wedge

The Counter Clockwise Wedge is used to create a counter clockwise rotation on a single tooth within an aligner.

UHFAPCC

Ref: 678-807




**Hu-Friedy**

### Clockwise Wedge

The Clockwise Wedge is used to create a clockwise rotation on a single tooth within an aligner.

UHFAPCW

Ref: 678-808



**Hu-Friedy**

### Eraser

The Eraser is designed to reverse or remove unwanted indents placed by other instruments in the Clear Collection.

UHFAPER

Ref: 678-811

# MAKE ME SMILE.





**GBT prevents demineralisation, white spots, caries, gingivitis, calculus formation and bad breath for orthodontic patients**

- Patients who wear orthodontic appliances are more likely to be subject to biofilm-related oral pathologies, such as caries or periodontal disease, as some areas in the mouth are more difficult to reach when brushing at home.
- Professional oral hygiene on a regular basis is therefore highly recommended to avoid complications.
- Not only does regular professional oral hygiene reduce the build-up of debris on and around brackets and attachments, but it also helps minimise wire friction.

**Minimally Invasive Prophylaxis on Orthodontic Patients  
Safe For All Orthodontic Appliances**



## METAL BRACKETS



## CLEAR ALIGNERS



**RELIANCE**  
**Ortho Flextech**

**Popular Choice**

- Flexible lingual retainer.
- Secure and stable long term passive retention.
- Adapts to lingual anatomy. Will take an active bend.
- Low profile for patient comfort.
- Stainless Steel offers increased tensile and torsional strength.
- Chairside application.



Stainless Steel Etched	UREL030
Stainless Steel Non Etched	UREL048
14kt White Gold Etched	UREL031
14kt White Gold Non Etched	UREL032
14kt White Gold Nickel Free	UREL033



Scan the QR Code: to watch the Reliance University clinical procedures.

**RELIANCE**  
**Bond-A-Braid Lingual Retainer Wire**



The dead soft wire easily adapts and is excellent for semi-permanent orthodontic splinting of the upper and lower incisors, diastema maintenance and fixed retention.

- Lingual retainer
- .027" x .011" ribbon arch
- 8-strand braided wire
- Each tube contains 10 (6") wires.

Pack 10 UREL034

**RELIANCE**  
**Retainium**



- Flexes Interproximally
- Lowest profile
- Flat .027" x .011" ribbon arch
- Adaptability
- MRI compatible
- Titanium Molybdenum Alloy
- Nickel Free.


Pack 10 UREL049

**RELIANCE**  
**Retainium XL**



- Flat .027" x .016" ribbon arch
- Low profile – comfortable and functional
- Increased strength
- MRI compatible
- Titanium Molybdenum Alloy
- Nickel Free.

Pack 10 UREL050



### Flattened Retainer Wire

Stainless Steel Flattened Twisted Retainer Wire.

Type	Wire Size	Code	Pack Size
Spool	Wire size .010 x .028	UMA0030	Each
Straight Lengths	Wire size .010 x .028	UMA0031	Pack of 10
Gold Plated	Wire size .010 x .028	UMA0032	Each
Tooth Coloured	Wire size .010 x .028	UMA0033	Each



everStick ORTHO is indicated for the aesthetic retention phase after active orthodontic treatment. The resin matrix holds the individual glass fibres in a bundle, which facilitates the handling and adaptation to the arch.

everStick ORTHO offers an aesthetic and cost-effective alternative for the retention phase after active orthodontic treatments.

- High aesthetics thanks to a good blending with the natural tooth structure.
- Easy to adapt directly on the teeth, eliminating the need for study models or second visits.
- Comfortable to wear and easy to clean.
- Completely metal-free, minimising the risk of allergies.
- Only everStick products have a unique, patented interpenetrating polymer network structure (IPN).

2 x 10cm Strips

UES08CM



**Kit contains:**

- 12 x Habit Control
- 12 x Cuspid Ramp
- 12 x Button Hook
- 12 x Molar Tube
- 12 x Bite Ramp 3 & 5mm
- 12 x Wire Bonder 3 & 5mm
- 2 x Brush handles.



**Reliamold Starter Kit**

Moulds are translucent and flexible to facilitate composite loading, curing and disengagement.

- A mould is filled with a light cure composite (GoTo®, T.K.O.®, or Flowtain™ recommended) then delivered to the tooth and cured.
- The desired auxiliary is delivered, formed, and functional immediately.

UREL023



**Button Attachment**

Moulds are translucent and flexible to facilitate composite loading, curing and disengagement.

25 Moulds per Package UREL024



**Reliamold Molar Pad Attachment**

Open the bite with posterior occlusal build-ups. T.K.O.® or Ultra Band-Lok® recommended.

25 Moulds per Package UREL025

**Attachments**



**Poly-Clear Buttons**

- Ultra Clear Plastic Buttons
- Super Low Profile
- Excellent for Aligner Treatment.
- Pack 10

Shape: Oval	UMA0010
Shape: Round	UMA0011
Shape: Trapezoidal	UMA0012



**Button Ceramic Lingual**

Pack 10.

UAST028



**GESTENCO Freddy Ceramic Button**

Pack 10.

UGSC001



**Bondable Button**

Pack 10.

Round	UAST026
Rectangular	UAST027



**GESTENCO Freddy Hook**


Pack 10.

UGSC002



**RELIANCE**  
**Perfect A Smile - Pontic Paint**

Quick and economical. Light Cure, apply smoothly to clear aligner sockets. Translucent to match corresponding enamel shade. High elasticity, bonds to any clear aligner.




4g Jar Bleach	UREL021
4g Jar A1 Shade	UREL012
4g Jar A2 Shade	UREL013
4g Jar A3.5 Dark Shade	UREL014
4g Jar B1 Shade	UREL015
4g Jar C1 Shade	UREL016
4g Jar C4 Shade	UREL017
4g Jar D3 Shade	UREL018
4g Jar Light Shade	UREL019
4g Jar Medium Shade	UREL020



**RELIANCE**  
**Perfect A Smile - Pontic Putty**

Finger tip pontic fabrication. Light cure composite paste with putty like consistency. Facilitates tooth shape formation. Eliminates pontics, wax & PVS.



5g Syringe	UREL022
------------	---------

Scan the QR Code: to watch the Reliance University clinical procedures.



### Startflow Flowable Composite (Starflow)

Startflow is a light-cure, highly filled flowable composite for anterior and posterior restorations, low wear applications, cavity liner and conservative restorations.

- Radiopaque
- Sub-micron particle size gives superior polishability and resistance to staining
- True Vita shade match
- Non-sticky
- Fluoride slow release and natural fluorescence.

A1	UEVI003
A2	UEVI001
A3	UEVI004
A3.5	UEVI002
A4	UEVI005
Replacement Tips	UEVI006



### VOCO Aligner Flow

Fluorescent composite for aligner attachments.

- Secure grinding and removal under UV-A light
- Precise application with perfect flow
- Light-curing through the transparent template
- Stable and abrasion resistant throughout the entire treatment period
- All Popular Vita Shades.

A1	2 x 2g	UALFSA1
A2	2 x 2g	UALFSA2
A3	2 x 2g	UALFSA3
44 Gauge Tips		UVDC044
46 Gauge Tips		UVDC046



### Flow Tain

- Needle Tip Dispenser
- Using tips provided, these materials are easily injected onto custom pads
- Low modulus of elasticity to resist breakage
- Highly filled to bridge gaps between custom pads and tooth surface
- Fluoresces under UV light

Flow Tain + 10 Tips	UREL042
Flow Tain LV + 10 Tips	UREL043
Flow Tain 19 gauge Tips x 50	UREL051





### GoTo Adhesive Paste in Push Syringe



Now Fluoresces Under Black Light for Complete Clean-up! GoTo is a light cure adhesive with superior handling properties.

It does not allow bracket flotation once placed on the tooth or during flash removal. A fraction of the retail cost of leading, competitive light cure pastes.

1 x 4gm GoTo Syringe	UREL036
10 x 0.4gm GoTo PLT Tips	UREL037
<b>New</b> GoTo Colour Changing Syringe 4gm	UREL128



### Assure Plus

All surface primer, eliminates hydrofluoric acid application. Enhances bond strength on all surfaces, tolerates moisture contamination. Primes and conditions any surface.



6ml Bottle	UREL035
------------	---------



### Bond Aligner Adhesive with Tips

Bonds Attachments To Clear Aligners

- Secure aesthetic bonds
- Bond to any clear aligner material
- No primers or conditioners
- Clear light cure paste
- High modulus of elasticity
- Conforms to aligner distortions
- Secure bonds throughout
- Bond ceramic, metal or plastic attachments.

1 x 1.4gm Syringe + 10 Tips	UREL038
-----------------------------	---------



### L.C.R. Paste in Push Syringe

L.C.R is a flowable, light cure orthodontic adhesive intended to be used within an orthodontic, dental or paediatric dental practice for the bonding of lingual retainers, creation of occlusal buildups and for retention of thermoplastic aligners.

5g L.C.R. Paste in Push Syringe Micro+ nano fill for flowable consistency with durability of highly filled composite. Smooth viscosity. Cures to a glassy finish. Patient comfort – reduced wear.

1 x 5gm Syringe	UREL039
1 x 5gm Syringe + 5 Tips	UREL040
11 x 0.45gm Tips	UREL041



# BioMin®

## Armour for Teeth

### Protect your teeth from white spot lesions during Orthodontic treatment with BioMin® F toothpaste



BioMinF Toothpaste 75ml	Pack 8	UBIOMINF
BioMinC Toothpaste 75ml	Pack 8	UBIOMINC

- Helps reduce the risk of white spot lesions during Orthodontic treatment
- Provides controlled release of Calcium, Phosphate and Fluoride for up to 12 hours after brushing
- Strengthens teeth and protects against decay
- Makes teeth feel cleaner and fresher
- Developed from research at Queen Mary's University of London
- Develops acid-resistant fluorapatite plugging of exposed dentine tubules.

[www.biomin.co.uk](http://www.biomin.co.uk)

@BioMinTech /BioMinTechnologiesLimited @biomintechnologieslimited



# Home Care Hygiene Range



## Clear Aligner Home Care Kit

This pack contains the essential items needed to help patients clean and care for their aligners as well as maintaining good Oral Hygiene.

UORC090

### Pack contents:

- Sonic cleaner
- Retainer Brite - aligner appliance cleanser (1 box 36 tablets)
- Interdental brushes (assorted pack)
- Aligner Chewies (help seat aligners correctly)
- Aligner remover
- Fresh & Go Mouthwash
- Aligner case
- Mirror
- All supplied in a convenient, large, heavy duty zip bag.



## CURAPROX Aligner Foam

Helps protect teeth underneath aligners, retainers and mouth guards. The foam restores protects and remineralises with natural minerals and forms a moisturising, protective film that feels silky smooth.

UCUR009 40ml



## Mini Aligner Home Care Kit

This pack contains the essential items needed to help patients clean and care for their aligners as well as maintaining good Oral Hygiene.

### Pack contents:

- 1 x Aligner Remover
- 2 x Retainer Brite Tablets
- 2 x Aligner Chewies
- 1 x Storage Case

UORC091



**Essix Retainer Brite**

Cleans and brightens retainers, night guards, mouth guards, splints and most other removable dental appliances.

- Helps remove plaque & tartar
- Kills germs & bacteria
- Keeps appliance minty-fresh.

Retainer Brite 36 - 1 month supply	URBT000
Retainer Brite 92 - 3 month supply	URBT001



**Retainer Boxes**

Assorted colours.

Small	Pack 24	UPRO002
* Small Colours: Red, Green, Pistachio, Blue, Fuchsia, Orange. 29x76x72mm		
Glitter Retainer Boxes	Pack 24	UPRO003
* Glitter Colours: Purple, Green, Yellow, Blue, Pink, Orange. 29x76x72mm		
Deep Pack	Pack 24	UPRO001
* Deep Colours: Red, Green, Pistachio, Blue, Fuchsia, Orange. 42x82x72mm		



**ReliaPry**

**Introducing ReliaPry Clear Aligner Removal Tool**

Designed to offer an effortless solution to the daily task of aligner removal.

With ReliaPry, your patients will have a reliable and gentle tool that allows removal of clear aligners easily and without stress.

Resealable and tamper-proof patient packs.

- Individually wrapped patient packs
- Each patient pack contains 1 ReliaPry tool
- Contains 100 patient packs - Assorted Colours.

Pack 100	UREL011
----------	---------

**Essix Art Appliance Remover Tool**



ART makes it easy for patients to remove their clear appliances! Removes clear appliances from the cheek-side and/or from the tongue side. Increases patient comfort and compliance. Latex free.

Pack 25	UORC100
---------	---------



**Retainer Box**

Bleaching/Retainer boxes in a variety of colours.

Shallow	Pack 10	UBRBSAC
* 88x73x30mm		
Deep	Pack 10	UBRBDSC
* 80x70x40mm		





## GC Tooth Mousse

To reduce white spot lesions and a topical application after IPR.

Intro Kit	Pack 5	UTMI000
Assorted	Pack 10	UTM0006
Mint	Pack 10	UTM0003
Strawberry	Pack 10	UTM0001



### Fluorid AF Daily Rinse

A daily mouth rinse that prevents dental caries. Suitable for children over 6 years of age.

Listed on the British National Formulary (BNF) as Sodium Fluoride mouth wash 0.05% (225ppm F) Mint flavour. Alcohol free.

Blue 400ml	Pack 4	UFAFDRB
Green 400ml	Pack 4	UFGFMRG

### Fresh & Go Mouthwash

Fresh & Go® Mouthwash is specifically formulated for orthodontic appliance wearers as it is alcohol free, contains fluoride and has a refreshing mild after taste.

- Alcohol free
- 0.05% Sodium Fluoride 225ppm F
- Made in the UK

500ml	Pack 12	UORC104
-------	---------	---------

## Clinisept+ Dental

Outstanding clinical results for mucosa.

### Clinisept

Clinisept+ Mouthwash is an oral cleanser and hygienic rinse that delivers professional levels of antimicrobial protection against the harmful bacteria that cause gum disease and tooth decay, but without stinging or causing irritation.

400ml	Pack 12	UCSPMW4
600ml	Pack 6	UCSPMW6

- Effective cleansing and antimicrobial protection (EN1040:2005, EN1276:2019, EN14476:2013 +A1:2015)
- Contains no chlorhexidine, alcohol, colouring or flavouring
- Non-irritant to teeth or gums
- Causes no staining
- Neutral pH 7.0
- Non-cytotoxic and hypo-allergenic

Clinisept+ Dental Mouthwash has an 18-month shelf life and should be used within 3 months of opening.



# BioCrystal Aesthetic Distalizer

Distalizer



- Simple Fabrication for Class II and III cases  
Bond the loose button to the first molar
- Measure the arm for bonding to the middle third of the canine allowing for full insertion into the first molar
- Bond a buccal tube to the opposing molar
- Run Class II or Class III Elastics to activate



## Bio-Crystal Distalizer Kit

- Stain Resistant Polycarbonate Aesthetics
- 4 sets of distalizers with tubes
- Fully Adjustable, clear ruler included
- Amazing clarity, function and very aesthetic
- 2 pairs of Distalizers Cost Effective
- Easy instructions included
- No Inventory Issues - One size fits all
- No Laboratory Fee.

Elastics for Activation.

Diameter	Force	Force Rating	Code
1/4"	6.5oz	F1	UREL046
3/16"	8.5oz	F2	UGHF476

Pack 4 UMA0034



**Essix  
C+ Retainer Sheets**

**Round 5"**

**Pack 100**

.040"	1mm	120mm	1C120-50	UORC064
-------	-----	-------	----------	---------

**Pack 50**

.040"	1mm	125mm	1C125-50	UORC065
-------	-----	-------	----------	---------

**Square 5" x 5"**

**Pack 50**

.040"	1mm	127mm	1C500-50	UORC063
-------	-----	-------	----------	---------

**Essix  
Plus Retainer Sheets**

**Round 5"**

**Pack 80**

.040"	120mm	1APL40C120	UORC081
-------	-------	------------	---------

.040"	125mm	1APL40C125	UORC082
-------	-------	------------	---------

**Square 5" x 5"**

**Pack 80**

.040"	5"	127mm	1APL40S5	UORC080
-------	----	-------	----------	---------

**Essix  
A+ Retainer Sheets**

**Round 5"**

**Pack 50**

.80"	2mm	120mm	11081	UORC041
------	-----	-------	-------	---------

.080"	2mm	125mm	11080	UORC042
-------	-----	-------	-------	---------

**Square 5" x 5"**

**Pack 50**

.030"	.75mm	10035-50	UORC030
-------	-------	----------	---------

.040"	1mm	10045-50	UORC034
-------	-----	----------	---------

.060"	1.5mm	10065	UORC038
-------	-------	-------	---------

.080"	2mm	10085	UORC040
-------	-----	-------	---------

**Essix  
Ace Retainer Sheets**

**Round 5"**

**Pack 50**

.030"	120mm	.75mm	1A320-50	UORC046
-------	-------	-------	----------	---------

.030"	125mm	.75mm	1A325-50	UORC047
-------	-------	-------	----------	---------

.040"	120mm	1mm	1A420-50	UORC050
-------	-------	-----	----------	---------

.040"	125mm	1mm	1A425-50	UORC049
-------	-------	-----	----------	---------

**Square 5" x 5"**

**Pack 50**

.030"	.75mm	1A350-50	UORC060
-------	-------	----------	---------

.040"	1mm	1A450-50	UORC048
-------	-----	----------	---------



**Perfection Plus Crown Scissors Curved**

Stainless steel Bee Bee Crown Scissors Curved, high quality, excellent functionality. Fully autoclavable.

UPPSCCC



### Whitening Gel - Day & Night Kits

**Boutique By Day** is a 6% hydrogen peroxide gel, containing potassium nitrate which is a proven desensitiser. Worn for only one hour per day, this gel is designed for patients who want a fast, comfortable and convenient solution to whiter teeth.

This gel releases large quantities of oxygen free radicals (which whiten the teeth) over a short period of time, therefore whitening the teeth, both quickly and safely. Fast results, which are visible in less than a week.

**Boutique By Night** is a carbamide peroxide gel available in a 10% or 16%, containing potassium nitrate, which is proven desensitiser. Worn for 4 hours, or overnight, this gel is designed for patients who would prefer to not have their custom trays in during the day and would prefer to whiten whilst they are asleep.

This gel releases oxygen free radicals (which whiten the teeth) over a longer period of time, due to its sustained release. This gradual approach will still give fantastic results, due the length of time the gel stays active.

<b>Day Kit 6% HP</b> 4 x 3ml	UBTK006	<b>Combo Kit 1x6% 3x16%</b> 4 x 3ml	UBCK134
<b>Night Kit 10% CP</b> 4 x 3ml	UBTK010	<b>Combo Kit 2x6% 2x16%</b> 4 x 3ml	UBCK224
<b>Day Kit 16% CP</b> 4 x 3ml	UBTK016	<b>Combo Kit 3x6% 1x16%</b> 4 x 3ml	UBCK314
<b>Night Combo Kit 1x10% CP &amp; 3x16% CP</b> 4 x 3ml	UBCK164	<b>Combo Kit 2x10% 1x16% 1x6%</b> 4 x 3ml	UBCK166
<b>Night Combo Kit 2x10% CP &amp; 2x16% CP</b> 4 x 3ml	UBCK264	<b>Combo Kit 2x16% 1x10% 1x6%</b> 4 x 3ml	UBCK106





### B1 Black Diamond 16% Kit

Kit Contains:

- Syringes 3ml x 8
- Whitening trays
- Boutique biomin toothpaste
- Boutique black charcoal floss
- B1 guarantee

UBTR010



### Whitening Gel Night Refill

10% or 16% carbamide peroxide gel with potassium nitrate. By Night is a 16% carbamide peroxide gel, containing potassium nitrate, which is proven desensitiser.

Pillow Box 10% CP	Pillow Box 16% CP
<b>1 x 3ml</b> UBPB101	<b>1 x 3ml</b> UBPB161
Pillow Box 10% CP	Pillow Box 16% CP
<b>2 x 3ml</b> UBPB102	<b>2 x 3ml</b> UBPB162
Pillow Box 10% CP	Pillow Box 16% CP
<b>3 x 3ml</b> UBPB103	<b>3 x 3ml</b> UBPB163



### Boutique Whitening Gel Night Hybrid Pro

Boutique whitening hybrid-pro 4.25% hydrogen peroxide 4.25% carbamide peroxide boutique hybrid-pro whitens your teeth overnight.

Pillow Box 4.25%	
<b>1 x 3ml</b>	UBHPPB1
Pillow Box 4.25%	
<b>2 x 3ml</b>	UBHPPB2
Pillow Box 4.25%	
<b>3 x 3ml</b>	UBHPPB3



### Whitening Gel Day Refill

Boutique – By Day is a 6% hydrogen peroxide gel, containing potassium nitrate which is a proven desensitiser. Worn for only one hour per day, this gel is designed for patients who want a fast, comfortable and convenient solution to whiter teeth.

Pillow Box 6% HP	
<b>1 x 3ml</b>	UBPB061
Pillow Box 6% HP	
<b>2 x 3ml</b>	UBPB062
Pillow Box 6% HP	
<b>3 x 3ml</b>	UBPB063

Whitening



360 whitening



360 whitening

### Whitening Gel - Day and Night Kits

The dual action of 360 Whitening gel lightens your patient's teeth from within while simultaneously removing surface stains.

Once applied to the teeth the bond gently breaks down creating free radicals that gently but effectively whiten the teeth quickly without sensitivity.



Day 6% HP Kit	<b>4 x 3ml</b>	U360064
Night 16% CP Kit	<b>4 x 3ml</b>	U360164
Night 10% CP Kit	<b>4 x 3ml</b>	U360104

### Whitening Refill Syringe

Day 6% HP Top Up	<b>1 x 3ml</b>	U360061	Day 6% HP Top Up	<b>3 x 3ml</b>	U360063
Night 10% CP Top Up	<b>1 x 3ml</b>	U360101	Night 10% CP Top Up	<b>3 x 3ml</b>	U360103
Night 16% CP Top Up	<b>1 x 3ml</b>	U360163	Night 16% CP Top Up	<b>3 x 3ml</b>	U360161



### BlancOne® HOME+ EVER Night-time whitening treatment for rejuvenation of teeth

Featuring  
**7%**  
carbamide  
peroxide\*



#### Indications

Orthodontic aligners  
Sensitive teeth  
Maintenance

More  
**BlancOne**  
choices  
online

#### Treatment cycle

10-15 applications - night time

Typical effectiveness

Stability

Regeneration

**BlancOne Home+ Ever Single Kit**

**Code: UBOHESK**

**BlancOne Home+ Ever Multi Kit**

**Code: UBOHEMK**

\*it can be used for extended periods.


Add an extra dimension to your Orthodontic treatment protocol with composite bonding

# Never puzzle with shade uncertainty again!



 Tokuyama **ESTELITE Σ QUICK**

Panoramic shades...  
...easily bridge the gap!




For simple restorations, rebuilding or layering, Estelite delivers every time!

E.g. 8 Shades replaced with A3  
in an anterior rebuild.



Scan the QR  
to view full  
case reports

# Available in 20 Vita Shades...



Estelite Sigma Quick has 20 Vita shades, varying in hue, chroma, opalescence and opacity.

Semi-opaque opalescent shades are designed for blocking out the darkness of the oral cavity in Class III & Class IV restorations.

The opaque shade is for masking slight stain or reconstructing highly opaque teeth (use as a thin layer at the lingual cavity wall of Class III & Class IV).

# ...but 3 Shades Can Do It All!

		B1	A1	B2	D2	A2	C1	C2	D4	A3	D3	B3	A3.5	B4	C3	A4	
Basic Shade (Class I, II, V)	A2	←————→															
	A3.5						←————→										
Opalescent (Class III, IV)	OA3			←————→													

Shade Matching Chart (Arrows indicate the range of tooth colour one resin shade can cover)  
 Syringe 3.8g PLT Pack 20 x 0.20g

# The 6 Shade Luxury Kit

contains: A1, A2, A3, A3.5, OA2, OPA2 & 25 x Microbrushes, 1 x 1ml Bond Force II, 1 x Shade Selector, 1 x Mixing Dish

A great way to get started!

		B1	A1	B2	D2	A2	C1	C2	D4	A3	D3	B3	A3.5	B4	C3	A4	
Basic Shade (Class I, II, V)	A1	←————→															
	A2	←————→															
	A3		←————→														
	A3.5			←————→													
Opalescent (Class III, IV)	OA2	←————→															
	OPA2		←————→														

Having 4 opaque shades (body), Estelite can be used as a multi-shade layering composite, ideal for Class IV.  
 Syringe Kit 6 x 3.8g. PLT Kit 40 x 0.2g.



To order or request a demo, contact us today.



**Trycare** [www.trycare.co.uk/dental](http://www.trycare.co.uk/dental)

Call: 01274 88 55 44 Email: [orthodontics@trycare.co.uk](mailto:orthodontics@trycare.co.uk)

Trycare Dental Rialto Buildings, Station Road, Clayton, Bradford, BD14 6JA, UK

■ ■ ■ **XXLITD146**

© 2025 Trycare Limited. All rights reserved. V2 01/25