

# 4 reasons to use Luxator® Short Periotome

## #1

### Ergonomic design

World acclaimed handle design. Minimizes fatigue. The 10mm shorter blade enables an optimal working position for smaller hands.



## #2

### Atraumatic extractions

The fine tapered blade compresses the alveolar bone and cuts the periodontal ligaments.

## #3

### Maximum tactility

The use of highest quality material gives the instruments outstanding tactility.

## #4

### Wide range of sizes and shapes

Five sizes and shapes enable users to find the ideal instrument for any extraction.



**S2S Periotome**  
2 mm  
Straight  
Short




**S3C Periotome**  
3 mm  
Curved  
Short

**S3S Periotome**  
3mm  
Straight  
Short

**S3CA Periotome**  
3 mm  
Contra  
Short

**S5S Periotome**  
5 mm  
Straight  
Short

## Technique



### Luxator® Short Periotome with Fixus1

Prod No	Product Description	Quantity
506358	S2S 2mm, straight blade, dark green	1 pc
506359	S3C 3mm, curved blade, black	1 pc
506360	S3S 3mm, straight blade, grey	1 pc
506361	S3CA 3mm, contra angle, dark blue	1 pc
506362	S5S 5mm, straight blade, tan	1 pc
<b>Kit</b>		
506333	Kit SK4 S3C, S2S, S3S, S3CA, L10	4 pcs