



Product Comparison

TABLETS	Team Vista	Mint-A-Clean	ICX/ Sterilex	BlueTab/ Sterilex	Citrisil™ / Citrisil™ Shock
Odorless	✓	✓	✓	✓	✓
Tasteless			✓	✓	✓
Non-Corrosive		✓	✓	✓	✓
Non-Toxic			✓	✓	✓
No need to empty bottles or purge lines				✓	✓
Included Shock Treatment					✓
Does not affect composite bonding					✓
EPA registered to disinfect dental output water ≤ 10 CFU/ml					✓
Option color indicator for treated water visual compliance					✓
Compatible with Antimicrobial Biofree™ or A-Dec tubing					✓
EPA registered to maintain Antimicrobial Biofree™ or A-Dec tubing					✓
STRAWS/ CARTRIDGES	Dentapure DP365B	Sterisil® Straw	Dentapure DP365M	Sterisil® Cartridge	
Odorless	✓	✓	✓	✓	
Tasteless		✓		✓	
EPA registered to disinfect dental output water ≤ 10 CFU/ml		✓		✓	
No need to empty and clean bottle or purge lines daily		✓		✓	
EPA registered shock treatment included		✓		✓	
Does not dissolve amalgam or cause possible discharged violations		✓		✓	
Removes source water total dissolved solids (TDS) to zero ppm		✓		✓	
Compatible with Antimicrobial Biofree™ or A-Dec tubing		✓		✓	
EPA registered to maintain Antimicrobial Biofree™ or A-Dec tubing		✓		✓	
SYSTEMS	O-So Pure	VistaClear	VistaPure	Sterisil® Systems	
Particulate Filtration	✓	✓	✓	✓	
Taste & Odor removal	✓	✓	✓	✓	
UV disinfection of broad range of bacteria & viruses				✓	
EPA registered residual disinfectant				✓	
Deionized autoclave water			✓	✓	
Drinking water				✓	
Automated system monitoring				✓	
EPA registered to disinfect dental output water ≤ 10 CFU/ml				✓	
Compatible with Antimicrobial Biofree™ or A-Dec tubing				✓	
EPA registered to maintain Antimicrobial Biofree™ or A-Dec tubing				✓	

Sterisil® , it doesn't get any easier, and it certainly doesn't get any better!

For more information about our products contact your **Sterisil® Sales Representative** or local dealer.

719 622 7200 or online at Sterisil.com